**CS Bezjehlový vstup M**

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**EN Needleless port M**

**Catalogue No.:** Medical device variants:  
 V696420-M Needleless port M  
 V696421-M Needleless port – red M  
**Functionality:** the manufacturer declares the functionality of the device when used in accordance with the instructions for use  
**Duration of use:** „Short term“ - continuous use for the time maximum 7 days or 100 connections  
**Design:** sterile, ethylene oxide sterilization  
**Target patient:** Patients undergoing infusion therapy.  
**Designated user:** Doctors, anesthetists, nurses, and other healthcare professionals in healthcare settings who perform infusion therapy.  
**Materials used:** PC, Silicone  
**Compatibility:** With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. The needleless port is intended for the application of liquid drugs or infusion solutions intended for injection or intravenous administration or for blood sampling using a syringe (without a needle). (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see [www.sukl.cz](http://www.sukl.cz) or [www.olecich.cz](http://www.olecich.cz)  
**Contraindication:** Do not use on any patient with a known allergic reaction to any of the product components.  
**Adverse effects:** No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.  
**Warning:** Do not use, if the sterile packing is opened or damaged. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.  
**Storage and transport:** The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!


**Intended use**

The needleless port is a closed safety system for use without a needle, ensuring microbiological safety, which is designed for repeated application as well as aspiration from intravenous accessories (syringes, infusion sets, cannulas, taps, catheters, extension hoses) with Luer type connection, in infusion therapy, injection therapy, sample taking, chemotherapy and it can also be used in magnetic resonance examinations. Red and blue color coding of the needleless ports makes it easy to identify the arterial or venous line. The replacement of the needle-free input must always be done after a maximum of 7 days or after 100 connections.

Priming volume: 0,09 ml  
 Flow Rate Averages: 1 psi: 550 ml/min = gravitational flow ; 3 psi: 1000 ml/min ; 5 psi: 1300 ml/min

**Instruction for use**

- Under aseptic conditions, remove the sterile needleless port from the packing.
- Prior to application, disinfect the silicone membrane surface by wiping with alcohol disinfection.
- Wait at least 30 seconds until the disinfection solution has dried out.
- By pushing on the Luer syringe or other Luer connector, carefully make the connection directly into the valve using a clockwise rotation hand to ensure a secure connection. **DO NOT** attempt to insert at an angle or try to pry the slot in the valve.
- Disconnecting the needleless port is carried out by turning anticlockwise.

**Replacement of the needleless port shall be performed always after 7 days or 100 connections at maximum.**  
 In order to prevent damage to the needleless port, use only the standard connection Luer which is in accordance with harmonized standards (non-standard syringes, needles, closures or connectors will cause damage to the needleless port). Do not angle.

**NOTE:** It is not necessary to use gloves while using the needleless port.

**Procedure for safe disposal of medical devices after use:**

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

**Used symbols on packaging**

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide and single sterile barrier system
	Pyrogen free		Keep away from sunlight		Latex free product		Unique Device Identifier
	Do not Resterilize		For single use only		Do not use if unit package is damaged		
	Medical Device		Date and Country of Manufacture		Does not contain PVC		

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