CS Bezjehlový vstup

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EN Needleless port

Catalogue No.: Medical device variants:
 V696420 Needleless port
 V696421 Needleless port – red

Product type: medical device

Design: sterile, ethylene oxide sterilization

Target patient: population of pediatric or adult patients

Designated user: medical personnel with appropriate education, knowledge and practical experience in medical devices.

Materials used: PC BPA-free / lipid resistant, Silicone

Compatibility: With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz

Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.

Adverse effects: No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.

Warning: Do not use, if the sterile packing is opened or damaged. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.

Storage: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.


Intended use

The needleless port is a closed safety system for use without a needle, ensuring microbiological safety, which is designed for repeated application as well as aspiration from intravenous accessories (syringes, infusion sets, cannulas, taps, catheters, extension hoses) with Luer type connection, in infusion therapy, injection therapy, sample taking, chemotherapy and it can also be used in magnetic resonance examinations.

Priming volume: 0,09 ml
 Flow Rate Averages: 1 psi: 550 ml/min = gravitational flow ; 3 psi: 1000 ml/min ; 5 psi: 1300 ml/min

Instruction for use

- Under aseptic conditions, remove the sterile needleless port from the packing.
- Prior to application, disinfect the silicone membrane surface by wiping with alcohol disinfection.
- Wait at least 30 seconds until the disinfection solution has dried out.
- Carefully connect (a common luer syringe or a luer lock connector) by turning clockwise to ensure a safe connection.
- Disconnecting the needleless port is carried out by turning anticlockwise.

Replacement of the needleless port shall be performed always after 7 days or 100 connections at maximum.
 In order to prevent damage to the needleless port, use only the standard connection Luer which is in accordance with harmonized standards (non-standard syringes, needles, closures or connectors will cause damage to the needleless port). Do not angle.

NOTE: It is not necessary to use gloves while using the needleless port.

Procedure for safe disposal of medical devices after use:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide
	Pyrogen free		Keep away from sunlight		Latex free product		
	Do not Resterilize		For single use only		Do not use if unit package is damaged		
	Medical Device		Date of Manufacture		Does not contain PVC		

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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