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Suction catheter M

Catalogue No.: Medical device variants:

V646738-M Suction catheter with CO, 1 BO, 2 mm (Ch 6) x 400 mm M
V646739-M Suction catheter with CO, 1 BO, 2,67 mm (Ch 8) x 400 mm M
V646740-M Suction catheter with CO, 1 BO, 3,33 mm (Ch 10) x 400 mm M

Catalogue No.: Medical device variants:

V646741-M Suction catheter with CO, 1 BO, 4 mm (Ch 12) x 400 mm M
V646742-M Suction catheter with CO, 1 BO, 4,67 mm (Ch 14) x 400 mm M
V646743-M Suction catheter with CO, 1 BO, 5,33 mm (Ch 16) x 400 mm M

Catalogue No.: Medical device variants:

V646744-M Suction catheter with CO, 1 BO, 6 mm (Ch 18) x 400 mm M
V646745-M Suction catheter with CO, 1 BO, 6,67 mm (Ch 20) x 400 mm M

Functionality: the manufacturer declares the functionality of the device when used in accordance with the instructions for use
Duration of use: „Short term“ - continuous use for the time maximum is 24 hours
Design: sterile, ethylene oxide sterilization
Target patient: Patients who need to suction secretions from the airways.
Designated user: Doctors, nurses and other healthcare professionals in healthcare settings who specialize in intensive care or emergency medicine and have a good knowledge of respiratory anatomy and physiology.
Materials used: PVC
Compatibility: With a suction tube that ends with a connector for vacuum control. The catheter are not intended for drug delivery.
Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects: No adverse effects are known if the catheter is used in accordance with the instructions specified in the Warning and Use sections.
Warning: Perform suction with short, intermittent vacuum for a maximum of 10-15 seconds.
Do not use, if the sterile packing is opened or damaged. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage and transport: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation.
Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action.
The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

Intended use

Suction catheters are single use medical devices intended for tracheal-bronchial suction from the lungs. The flexible tube is designed to be inserted into the airways for the purpose of removing foreign material using suction.

Instruction for use

1. Check the type labelling of the catheter, usable life on the package and the package for damages.
2. Take the catheter out of the packaging under aseptic condition.
3. Insert the catheter into patient's airways.
4. Connect the catheter ending to the suction device.
5. Then follow the instructions for use of the suction device for which the catheter is intended.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide and single sterile barrier system
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Medical Device
	Date and Country of Manufacture		Unique Device Identifier				

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