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EN Catheter for neonates of low birth weight CN-01 M

Catalogue No.: V646959-M
Functionality: the manufacturer declares the functionality of the device when used in accordance with the instructions for use
Duration of use: „Short term“ - continuous use for the time maximum is 24 hours
Design: sterile, ethylene oxide sterilization
Target patient: Newborns with low birth weight (NLBW), children with a birth weight below 2500 g.
Designated user: Neonatologists in medical facilities who have relevant education and experience in the care of newborns.
Materials used: PVC
Compatibility: With Luer single-use syringe, Luer Lock single-use syringe and Infusion set (see Manufacturer's note). (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz
Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects: No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning: The catheter shall not be used to evacuate the air passages.
 Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date.
 Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage and transport: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation.
 Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action.
 The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

Intended use

Universal tube designed for neonates with low birth weight and allows use both for enteral nutrition and for infusions and application of medicine into the umbilical vein. It is a single-use medical device.

Instruction for use

1. Enteral nutrition fed to infants:
 Take the catheter out of the packing under aseptic condition and introduce it the rounded tip first to the infant's stomach through its mouth or nose. The depth of catheter penetration can be checked using the graduation marks of 5, 10, 15 and 20 cm as measured from the rounded tip. Before the nutritious solution can be administered, make sure the catheter did not enter the lower air passages. Then attach the catheter via the connecting cone to a syringe or an infusion set and start delivering the solution.

2. Medication to the umbilical vein
 Short-time infusions to the umbilical vein:
 Take the catheter out of the packing; open the cone closure and attach the cone to a syringe with infusion solution; fill the catheter with the solution. The rounded tip first introduce the catheter to the patient's umbilical vein. The depth of catheter penetration can be checked using the graduation marks of 5, 6 and 7 cm as measured from the rounded tip. Apply a short suction or relief action to make sure the catheter end is being filled with blood. Then re-attach the catheter to the syringe, the infusion set, etc., as the job at hand requires.

Manufacturer's note: the manufacturer warns the user of the existing possibility of incorrect connection. The connection connector used on the catheter is a separately developed alternative connector. Compatibility with syringe and infusion set has been verified by the manufacturer in physical and functional tests.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:
 Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Umbrella	Keep away from Rain	LOT	Batch Number	Factory	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide and single sterile barrier system
Pyrogen	Pyrogen free	Sun	Keep away from sunlight	Latex	Latex free product	PHT/DEHP	Does not contain DEHP phthalate
Resterilize	Do not Resterilize	Single use	For single use only	Damaged	Do not use if unit package is damaged	MD	Medical Device
Date	Date and Country of Manufacture	UDI	Unique Device Identifier				

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