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## Neonatal nutrition catheter CV-01 M

<b>Catalogue No.:</b>	V646957-M
<b>Functionality:</b>	the manufacturer declares the functionality of the device when used in accordance with the instructions for use
<b>Duration of use:</b>	„Short term“ - continuous use for the time maximum is 24 hours
<b>Design:</b>	sterile, ethylene oxide sterilization
<b>Target patient:</b>	A newborn (neonatalogical patient) who has a feeding problem.
<b>Designated user:</b>	Neonatologists, nurses and other healthcare workers in healthcare facilities who have relevant education and experience in caring for newborns.
<b>Materials used:</b>	PVC, PP
<b>Compatibility:</b>	With Luer single-use syringe, Luer Lock single-use syringe and infusion set (see Manufacturer's Note). CV-01 is not intended for drug delivery
<b>Contraindication:</b>	Do not use on any patient with a known allergic reaction to any of the product components.
<b>Adverse effects:</b>	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
<b>Warning:</b>	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
<b>Storage and transport:</b>	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

### Intended use

A single-use medical device intended for enteral feeding of neonates.

### Instruction for use

Take the food-delivery catheter out of the packing under aseptic condition and introduce it to the infant's pharynx the rounded tip with the side holes first. The depth of penetration can be checked using the graduation marks of 5, 10, 15 and 20 cm. Open the closure of the negative connecting cone and attach the catheter to the connecting cone of the infusion set; start delivering the liquid food.

**Manufacturer's Note:** the manufacturer warns the user of the existing possibility of incorrect connection. The connection connector used on the catheter is a separately developed alternative connector. Compatibility with syringe and infusion set has been verified by the manufacturer in physical and functional tests.

#### PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

### Used symbols on packaging

REF	Catalogue Number	Rain	Keep away from Rain	LOT	Batch Number	Factory	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide and single sterile barrier system
Pyrogen	Pyrogen free	Sun	Keep away from sunlight	Latex	Latex free product	DEHP	Does not contain DEHP phthalate
Sterilize	Do not Resterilize	Single use	For single use only	Damaged	Do not use if unit package is damaged	MD	Medical Device
Date and Country	Date and Country of Manufacture	UDI	Unique Device Identifier				

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