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**EN Neonatal nutrition catheter CV-01**

**Catalogue No.:** V646957-ND  
**Medical device variants:** Neonatal nutrition catheter CV-01  
**Product type:** medical device  
**Design:** sterile, ethylene oxide sterilization  
**Target patient:** newborn (neonatal patient)  
**Designated user:** medical personnel with appropriate education, knowledge and practical experience in medical devices.  
**Materials used:** PVC, PP  
**Compatibility:** Infusion sets and syringes must be equipped with Luer connection cone in conformity with standard ČSN EN ISO 80369-7. CV-01 is not intended for drug delivery  
**Contraindication:** Do not use on any patient with a known allergic reaction to any of the product components.  
**Adverse effects:** No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.  
**Warning:** Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.  
**Storage:** The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

**Intended use**

Enteral nutrition fed to infants.

**Instruction for use**

Take the food-delivery catheter out of the packing and introduce it to the infant's pharynx the rounded tip with the side holes first. The depth of penetration can be checked using the graduation marks of 5, 10, 15 and 20 cm. Open the closure of the negative connecting cone and attach the catheter to the connecting cone of the infusion set; start delivering the liquid food.

**PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:**

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

**Used symbols on packaging**

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide
Pyrogen free		Keep away from sunlight	Latex free product	Does not contain DEHP phthalate	
Do not Resterilize		For single use only	Do not use if unit package is damaged	MD	Medical Device
Date of Manufacture					

The instruction for use is also available in electronic form at [www.gama.cz](http://www.gama.cz). Free copies of other copies are available on request via e-mail: [sales@gama.cz](mailto:sales@gama.cz).

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