Cévka pro výživu novorozence CV-01

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EN Neonatal nutrition catheter CV-01

Catalogue No.: V646957-ND
Medical device variants: Neonatal nutrition catheter CV-01
Product type: medical device
Design: sterile, ethylene oxide sterilization
Target patient: newborn (neonatal patient)
Designated user: medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used: PVC, PP
Compatibility: With disposable syringe or infusion set (see Manufacturer's Note). CV-01 is not intended for drug delivery
Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects: No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning: Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

Enteral nutrition fed to infants.

Instruction for use

Take the food-delivery catheter out of the packing under aseptic condition and introduce it to the infant's pharynx the rounded tip with the side holes first. The depth of penetration can be checked using the graduation marks of 5, 10, 15 and 20 cm. Open the closure of the negative connecting cone and attach the catheter to the connecting cone of the infusion set; start delivering the liquid food.

Manufacturer's Note: the manufacturer warns the user of the existing possibility of incorrect connection. The connection connector used on the catheter is a separately developed alternative connector. Compatibility with syringe and infusion set has been verified by the manufacturer in physical and functional tests.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

Catalogue Number	Keep away from Rain	Batch Number	Manufacturer
Use by / Expiry date	Temperature Limitation	Read instructions for use	Sterilised by Ethylene Oxide
Pyrogen free	Keep away from sunlight	Latex free product	Does not contain DEHP phthalate
Do not Resterilize	For single use only	Do not use if unit package is damaged	Medical Device
Date of Manufacture			

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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