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Neonatal umbilical catheter CP-01

Catalogue No.:	V646958-ND
Product type:	medical device
Design:	sterile, ethylene oxide sterilization
Target patient:	neonate
Designated user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used:	PVC, PP
Compatibility:	With disposable syringe or infusion set (see Manufacturer's note). (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects:	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning:	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

- Medication to the umbilical vein (antibiotics, pressure-boosting drugs, respiratory aids, etc.)
- Short-time infusion to the umbilical vein applied when resuscitating the infant
- Evacuation of air passages (incl. the lower air passages)

Instruction for use

to items 1 & 2 Medication to the umbilical vein

Short-time infusion to the umbilical vein applied when resuscitating the infant:

Take the umbilical catheter out of the packing under aseptic condition; open the cone closure and connect the cone to a syringe with infusion solution; fill the catheter with the solution. The rounded tip first introduce the catheter to the patient's umbilical vein. The depth of catheter penetration can be checked using the graduation marks of 5, 6 and 7 cm as measured from the rounded tip. Apply a short suction or relief action to make sure the catheter end is being filled with blood. Then re-attach the catheter to the syringe, the infusion set, etc., as the job at hand requires.

to item 3 Evacuation of air passages (incl. the lower air passages):

Take the umbilical catheter out of the packing; open the cone closure and introduce the catheter to the infant's air passage the rounded tip first. Then attach it to the suction-generating device.

Manufacturer's note: the manufacturer warns the user of the existing possibility of incorrect connection. The connection connector used on the catheter is a separately developed alternative connector. Compatibility with syringe and infusion set has been verified by the manufacturer in physical and functional tests.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide
Pyrogen free		Keep away from sunlight	Latex free product	Does not contain DEHP phthalate	
Do not Resterilize		For single use only	Do not use if unit package is damaged	MD	Medical Device
Date of Manufacture					

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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