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## Rectal catheter M

<b>Functionality:</b>	the manufacturer declares the functionality of the device when used in accordance with the instructions for use
<b>Duration of use:</b>	„Short term“ - continuous use for the time maximum is 24 hours
<b>Design:</b>	sterile, ethylene oxide sterilization
<b>Target patient:</b>	Patients who need to inject drugs into the anus, perform an enema or examine the anus.
<b>Designated user:</b>	Medical personnel such as doctors, nurses or midwives who have sufficient knowledge and experience with this procedure to avoid harm to the patient's health.
<b>Materials used:</b>	PVC
<b>Compatibility:</b>	With a Janette lavage syringe and an irrigator for introducing enemas and rinsing. With rectal drugs and solutions. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see <a href="http://www.sukl.cz">www.sukl.cz</a> or <a href="http://www.olecich.cz">www.olecich.cz</a> .
<b>Contraindication:</b>	Do not use on any patient with a known allergic reaction to any of the product components.
<b>Adverse effects:</b>	No adverse effects are known if the catheter is used in accordance with the instructions specified in the Warning and Use sections.
<b>Warning:</b>	Do not use, if the sterile packing is opened or damaged. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
<b>Storage and transport:</b>	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

### Intended use

A single use medical device intended for rectal examination, medication administration and rectal drainage.

### Instruction for use

1. Check the type labelling of the catheter, usable life on the package and the package for damages.
2. Take the catheter out of the packaging under aseptic condition.
3. Use the catheter for the intended purpose.

### PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

### Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use	Sterilised by Ethylene Oxide and single sterile barrier system	
Pyrogen free		Keep away from sunlight	Latex free product	Does not contain DEHP phthalate	
Do not Resterilize		For single use only	Do not use if unit package is damaged	Medical Device	
Date and Country of Manufacture		UDI	Unique Device Identifier		

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