CS Combi zátky LL

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EN LL combo stoppers

Catalogue No.: Z676004, Z676005, Z676006
Medical device variants: LL combo stopper – red NO PVC, LL combo stopper – blue NO PVC, LL combo stopper – white NO PVC
Product type: medical device
Design: sterile, ethylene oxide sterilization
Target patient: population of pediatric or adult patients
Designated user: medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used: PP
Compatibility: With medical devices equipped with standard Luer Lock connectors in conformity with standard ČSN EN ISO 80369-7. LL combo stopper is not intended for drug delivery.
Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects: No adverse effects are known if the combo stopper is used in accordance with the instructions specified in the Warning and Use sections.
Warning: Do not use, if the sterile packing is opened or damaged. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

The LL combo stopper is designed for safety closures, where a Luer cap is required. The combination of MLL and FLL connecting beams secures a safe and versatile closure for syringes, intravenous catheters, three-way valves, ramps and other infusion systems. Colour coding of the combo stoppers makes it easy to identify individual infusion lines or therapeutic applications. The LL combo stoppers are designed for both pressure and gravity systems up to a pressure of 2 bar and are compatible with Luer-Lock cones, which meet the requirements of ČSN EN ISO 80369-7.

Instruction for use

- Under aseptic conditions remove the stopper from the packaging.
- Attach the stopper with a clockwise twisting motion so as to ensure a safe connection.
- Detach the stopper with a counter clockwise twisting motion.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use		Sterilised by Ethylene Oxide
Pyrogen free		Keep away from sunlight	Latex free product		Does not contain DEHP phthalate
Do not Resterilize		For single use only	Do not use if unit package is damaged		Medical Device
Date of Manufacture					

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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