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## Bloodline tubing system DIS 06-16® UNIV

<b>Catalogue No.:</b>	V606016-02
<b>Product type:</b>	medical device
<b>Design:</b>	sterile, ethylene oxide sterilization
<b>Target patient:</b>	population of pediatric or adult patients
<b>Designated user:</b>	medical personnel with appropriate education, knowledge and practical experience in medical devices.
<b>Materials used:</b>	PVC, PE, PP, rubber
<b>Compatibility:</b>	With end pieces of haemodialysers in conformity with standard ČSN EN ISO 8637. With end pieces of vascular access, stress gauges, heparin tubes having threads in conformity with standard ČSN EN ISO 80369-7. The bloodline tubing system DIS 06-16 UNIV is used to transfer blood, substitution, dialysis or saline solution or heparin. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz.
<b>Contraindication:</b>	Do not use on any patient with a known allergic reaction to any of the product components.
<b>Adverse effects:</b>	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
<b>Warning:</b>	Before using this product, you should carefully read the instructions for use and operation of the dialysis monitor. Both before and during treatment, you should make sure that all connections are tight and the set is properly established in the dialysis monitor. Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
<b>Storage:</b>	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

### Intended use

Dialysis sets are intended for single use with appropriate dialysis monitor for extracorporeal blood lavage:  
GAMBRO AK 100, 200, 200S, 95, 96, 100 ULTRA, 200 ULTRA, 200 S ULTRA  
FRESENIUS 4008  
Nikkiso DBB-05, DBB-07  
B-BRAUN Dialog, Dialog+

Dialysis sets can be used according to physician's directions only.

Before the use of the medical device read carefully instructions for use and follow the directions and notes here given, follow the directions of dialysis monitor and dialysis machine manufacturer, follow the instructions for use for other medical devices used in dialysis treatment respectively. Do not use the medical device in case the packing is damaged or when the covering caps have fallen down. Before and during dialysis treatment make sure that all connections are tight and set is correctly inserted into various parts of dialysis monitor. The air penetrating into extracorporeal blood circulation can cause air embolism.

**Insert pump segment of dialysis sets into blood pump immediately before flushing and dialysis process.**

**Not doing so may cause mechanical deformation of pump segment and after activation of blood pump may cause its perforation.**

### GENERAL INFORMATION

Type designation	Arterial set length	Venous set length	Arterial set volume	Venous set volume	Pump segment*)
DIS 06-16® UNIV	4,01 m	3,08 m	110 ± 3 ml	80 ± 3 ml	8/12 x 340 mm, tl. 2 mm

\*) Meaning of numbers used in pump segment: 8.0 mm inner diameter, 12.0 mm outer diameter x segment length, segment wall thickness.

Arterial part of set is in red colour, venous part is in blue colour. Maximum pressure 500 mmHg.

Flexible PVC of 61 ± 3 ShA hardness is used as material for all pump segments. Other parts which are in direct contact with blood are made of PVC with various hardness. Parts which are not in direct contact with blood are made of PE, PP, etc.

### Instruction for use

- Take out the arterial and venous part of the set and collecting bag from the unitare package immediately before the use. Follow aseptic way of use when handling the set.
- Connect set, using the connectors at the end of the set, to appropriate haemodialyzer. Screw the red connector to red (arterial) haemodialyzer output, screw the blue connector to blue (venous) haemodialyzer output. Follow all instructions given in monitor instructions for use when connecting the set.
- Insert pump segment into blood pump of dialysis monitor, follow instructions for use for the appropriate type of dialysis monitor and check correct segment insertion and pump adjustment.
- Flush haemodialyzer with connected dialysis set with saline solution and follow operating instructions for the given monitor type. Saline solution which is pushed out by blood can be collected into collecting bag.
- Use only approved needles or catheters with appropriate end and diameter when connecting the set to patient.
- Adjust the levels at 2/3 – 3/4 of the volume of dripping chambers at arterial and venous branch. Adjust the level using the tube and sterile syringe. The tube must be permanently closed during the haemodialysis treatment; it can be opened during level adjustment only. Venous chamber level must be always higher than detector active zone level.
- Connect tube for heparin administration to heparin pump. Check the tightness and security of connection. Anticoagulant treatment can be administered according to physician's directions. Regularly check heparin pump flow during haemodialysis in order to guarantee the dosage accuracy.
- Insert the pad (underpressure detector) if provided into its holder; press the switch first and insert the pad into holder clamp following the manufacturer's instructions for use.
- Measure arterial and venous blood pressure using appropriate connectors only, provided with separating filters.
- Injection parts at both arterial and venous parts of the set are intended for drugs administration. Use max.21G sterile needle for injection part membrane needling.
- Reverse flush is performed by routine procedure according to physician's directions; follow instructions given by monitor manufacturer.

### NOTE:

The manufacturer is not responsible for use of dialysis sets with inappropriate dialysis monitors, for use which is not according to instructions for use, for failures caused by inappropriate storage and handling. The manufacturer is also not responsible for user neglect when checking the condition of dialysis set before the use.

### Procedure for safe disposal of medical devices after use:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

### INFORMATION FOR THE RISK GROUP OF PATIENTS

Should treatment using this medicinal product involve treatment of children, pregnant or breast-feeding women, please become familiar with the following possible residual risks which may arise if using the product.

### Residual risks with the application of this medicinal product

The product contains phthalates. If phthalates are inhaled, consumed or penetrate the skin, they may cause cancer or increase its incidence. They may also cause congenital genetic damage or increase its incidence and may cause or increase the incidence of non-congenital negative effects on offspring or exacerbate the male or female reproductive functions or abilities.

### Explanatory notes to the symbols on the unit container

PHT – medicinal product containing phthalates

DEHP – di-2-ethylhexyl phthalate is contained in the medicinal product in a quantity of more than 0.1 hm. %.

### Information about the risk of the repeated use of a single-use medicinal product

The manufacturer does not guarantee the safety of the product if used repeatedly.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

### Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use		Sterilised by Ethylene Oxide
Pyrogen free		Keep away from sunlight	Attention		Contains phthalate DEHP
Do not Resterilize		For single use only	Do not use if unit package is damaged		Medical Device
Date of Manufacture					

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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