CS

GAMAPLUS 3,0 x 150 LCH M / GAMAPLUS 4,5 x 20 LCH M

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GAMAPLUS 3,0 x 150 LCH M / GAMAPLUS 4,5 x 20 LCH M

| | |
|-------------------------------|---|
| Catalogue No.: | V606372-M / V606372-01M |
| Functionality: | the manufacturer declares the functionality of the device when used in accordance with the instructions for use |
| Duration of use: | „Short term“- continuous use for the time maximum is 24 hours |
| Design: | sterile, ethylene oxide sterilization |
| Target patient: | Patients undergoing intravenous or intra-arterial infusion (transfusion) therapy where an extension or flexible connection between the infusion (transfusion) system and the intravenous/intra-arterial catheter is required. |
| Designated user: | Doctors, anesthesiologists, nurses and other healthcare workers in healthcare facilities who perform intravenous or intra-arterial infusion (transfusion) therapy. |
| Materials used: | PVC, PE, ABS |
| Compatibility: | With universal gradual (step) connectors. With Female Luer Lock connector which is in accordance with the standard ČSN EN ISO 80369-7. GAMAPLUS is designed to transfer liquid or i.v. administration or transfer of blood. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz |
| Contraindication: | Do not use on any patient with a known allergic reaction to any of the product components. |
| Adverse effects: | No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections. |
| Warning: | Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient. |
| Storage and transport: | The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days! |

Intended use

A single-use medical device designed to connect two different transfer sets which are ended with the negative cone and the gradual connecting cone.

Instruction for use

- Under aseptic conditions, remove the item from the packing.
- Remove the positive cone cover.
- Connect the positive cone to the negative cone of the transfer set.
- Connect the white terminal to the gradual connecting cone of the transfer set.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

| | | | | | | | |
|--------------------|----------------------|---------------------------------|---------------------------------|-------------------|---------------------------------------|--------------|--|
| REF | Catalogue Number | Rain | Keep away from Rain | LOT | Batch Number | Manufacturer | Manufacturer |
| Hourglass | Use by / Expiry date | Temperature | Temperature Limitation | Read instructions | Read instructions for use | STERILE EO | Sterilised by Ethylene Oxide and single sterile barrier system |
| Pyrogen free | Pyrogen free | Sunlight | Keep away from sunlight | Latex free | Latex free product | DEHP | Does not contain DEHP phthalate |
| Do not Resterilize | Do not Resterilize | Single use | For single use only | Damaged | Do not use if unit package is damaged | P | Pressure |
| MD | Medical Device | Date and Country of Manufacture | Date and Country of Manufacture | VOL | Volume of tubing | UDI | Unique Device Identifier |

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