CS GAMAPLUS M

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EN **GAMAPLUS M**

Catalogue No.: V606301-M V606301-01M V606303-M V606304-M V606304-01M	Medical device variants: GAMAPLUS 1,8 x 450 LL M GAMAPLUS 1,8 x 450 LLK M GAMAPLUS 1,8 x 1500 LL M GAMAPLUS 1,8 x 1800 LL M GAMAPLUS 1,8 x 1800 LLK M	Catalogue No.: V606306-M V606307-M V606321-M V606342-M	Medical device variants: GAMAPLUS 1,8 x 450 UNIV M GAMAPLUS 1,8 x 1800 UNIV M GAMAPLUS 3,0 x 150 LL M GAMAPLUS 3,0 x 300 L ML	Catalogue No.: V606347-M V606349-M V606371-M	Medical device variants: GAMAPLUS 3,0 x 2500 LL M GAMAPLUS 3,0 x 3500 LL M GAMAPLUS 3,0 x 150 YLL M
Functionality: the manufacturer declares the functionality of the device when used in accordance with the instructions for use					
Duration of use: „Short term“. If blood, blood derivatives, fat emulsions are administered or if the line had to be disconnected, continuous use is for a maximum of 24 hours, in case of administration of the same type of drug, use is for a maximum of 72 hours, in a closed system terminated with a needle-free inlet the interval can be extended up to 7 days. Immediate replacement is necessary in case of contamination of the line and in case of incompatibility. Replacement according to national guidelines and/or procedures of the healthcare facility. Depends on the intended treatment according to the summary of product characteristics of the drug/solution. sterile, ethylene oxide sterilization					
Design: Patients undergoing intravenous or intra-arterial infusion (transfusion) therapy where an extension or flexible connection between the infusion (transfusion) system and the intravenous/intra-arterial catheter is required.					
Target patient: Doctors, anesthesiologists, nurses and other healthcare workers in healthcare facilities who perform intravenous or intra-arterial infusion (transfusion) therapy.					
Designated user: PVC, PE					
Materials used: With medical devices equipped with standard Luer connectors in conformity with standards ČSN EN ISO 80369-7. GAMAPLUS is designed to transfer liquid or i.v. administration or transfer of blood. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz					
Compatibility: Do not use on any patient with a known allergic reaction to any of the product components.					
Contraindication: No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.					
Adverse effects: Use up to a pressure of max. 2 bar. Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date.					
Warning: Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.					
Storage and transport: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!					

Intended use

A single use medical device designed to extend or connect one or two infusion or transfusion sets to an intravenous needle..

Instruction for use

1. Take GAMAPLUS out of the packing under aseptic conditions.
2. Remove the covers.
3. Attach the negative (female) cone to the positive (male) cone of the set and fill the hose with liquid.
4. Connect the positive cone to the intravenous needle or a catheter.

Procedure for safe disposal of medical devices after use:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Rain	Keep away from Rain	LOT	Batch Number	Manufacturer	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide and single sterile barrier system
Pyrogen	Pyrogen free	Sun	Keep away from sunlight	Latex	Latex free product	DEHP	Does not contain DEHP phthalate
Sterilize	Do not Resterilize	Single use	For single use only	Damaged	Do not use if unit package is damaged	Pressure	Pressure
VOL	Volume of tubing	MD	Medical Device	Date	Date and Country of Manufacture	UDI	Unique Device Identifier

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