CS HDF on-line set

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Catalogue No.:	V686435-ND
Product type:	medical device
Design:	sterile, ethylene oxide sterilization
Target patient:	population of pediatric or adult patients
Designated user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used:	PVC, PE, PP, PC, acrylic, silicone rubber
Compatibility:	With dialysis monitor DIALOG PLUS. With dialysis sets in conformity with standard ČSN EN ISO 8638. HDF on-line set is not intended for drug delivery.
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects:	In case of improper handling there is danger of a significant increase in occurrence of bacteria or quantity of pyrogene in the set. This results in threat of adverse reactions in the patient body. In this case, it is necessary to immediately stop the inlet blood flow and provide adequate medical care.
Warning:	Before using this product, you should carefully read the instructions for use and operation of the dialysis monitor. Both before and during treatment, you should make sure that all connections are tight and the set is properly established in the dialysis monitor. Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.


Intended use

Disposable medical device intended for application of substitution solution for the process of on-line hemo-diafiltration which is to be prepared in the appliance during the procedure performed on dialysis monitors of Dialog Plus made by B-Braun.

Instruction for use

You should establish the HDF on-line set on the appliance during its preparation only after installing the blood sets, but always before starting the automatic test of the appliance.

- Take the set out of the package.
 - Open the cover of the substitution (venous) pump. Insert the pump segment into the pump and then close the cover.
 - Make sure that the set is established in the pump in the right direction (the pump should rotate in the clockwise direction during operation).
 - Open the output of the substitution port (right, white) and connect it to the input part of the HDF on-line set.
 - Then you should connect the output part of the HDF on-line set to the venous bag of the blood set while in the process of post dilution HDF. While the pre-dilution HDF is in progress, you should connect it to the connector of the arterial blood set before the dialyzer.
 - Then you should connect the patient's end of the arterial blood set to the Luer-lock connector of this HDF on-line set between the substitution port of the appliance and the substitution pump.
 - Patient's end of the venous blood set should be connected to the left (blue) substitution port of the appliance.
- During all the time of the above-described manipulation, you should take care not to contaminate any connector of the sets or substitution ports of the appliance. To start flushing the entire extra-corporeal circuit and subsequent connection of the patient, including start and end of the HDF, you should comply with the instructions for operating the appliance.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
	Use by / Expiry date	Temperature Limitation		Read instructions for use	Sterilised by Ethylene Oxide
	Pyrogen free	Keep away from sunlight		Latex free product	Does not contain DEHP phthalate
	Do not Resterilize			Do not use if unit package is damaged	Medical Device
	Date of Manufacture				

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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