CS Infuzní souprava IS-121 M, IS-123 M

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EN Infusion set IS-121 M, IS-123 M

Catalogue No.: V606121-M, V606123-M
Functionality: the manufacturer declares the functionality of the device when used in accordance with the instructions for use
Duration of use: „Short term“ - continuous use for the time maximum is 24 hours. Is dependent on the intended therapy according to the SPC of the drug/solution. Change the set according to national guidelines and/or hospital protocols.
Design: sterile, ethylene oxide sterilization
Target patient: Patients undergoing intravenous or intra-arterial infusion therapy.
Designated user: Physicians, anesthetists, nurses, and other healthcare professionals in healthcare settings who administer intravenous or intra-arterial infusion therapy.
Materials used: IS-121 M = PVC, PE, ABS, PP, PS, NYLON, Membrane VERSAPORE RC, POLYISOPRENE; IS-123 M = PVC, PE, ABS, PP, PS, NYLON, membrane VERSAPORE RC, POLYISOPRENE, SILICONE
Compatibility: With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With plastic infusion bags, plastic or glass bottles. With infusion pharmaceuticals without any special requirements for filtration. (In)compatibility of the drugs with the materials used check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz
Incompatibility: With infusion pharmaceuticals incompatible with PVC (cytostatics) and with infusion pharmaceuticals sensitive to light.
Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects: No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning: Do not use if package is open, damaged, guards are dropped or have expired. This may compromise the integrity or sterility of the device. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It can lead to contamination or impaired functionality, which can cause injury, illness, or death to the patient. Connect the infusion bottle to the infusion set just before administration. During the procedure, check the patient's condition and the functionality of the set. Use the set only in conjunction with the prescribed pump type!
Storage and transport: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

Intended use
 Infusion sets for positive pressure infusions are disposable medical devices intended for the therapeutic administration of infusion solutions from bottles or bags using infusion pumps, with each type of set intended for a specific type of infusion pump:
 IS-121 M infusion set: Star-flow 581 infusion pump / IS-123 M infusion set: Star-flow 591 infusion pump.
When using the sets see to it that the attached pumps are of the stipulated type.

Instruction for use
 1. Under aseptic conditions, remove the set from the packing.
 2. Close the plug in the aeration hole on the body of the piercing point.
 3. Close the control clamp on the hose downstream of the drip chamber.
 4. Detach the cap from the piercing point and remove the cone cover.
 5. Insert the piercing point through the stopper to a bottle or a bag. Hang the bottle or the bag on a rack.
 6. Press the drip chamber so that it fills up to one half, open the control clamp, open the plug in the aeration hole and drive air out of the hose.
 Note: In case the set is to transport saline from the plastic bag, keep the aeration hole closed.
 7. Close the control clamp.
 8. Then follow the instructions on use of the pump to be installed jointly with the set.
NOTE: With the pump operational, the control clamp must be fully open. With the dosing process finished, the control clamp is to be closed before the set may be removed from the device.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:
 Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.
NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

| | | | | | | | |
|--|------------------------------------|--|---------------------------------|--|---------------------------------------|--|--|
| | Catalogue Number | | Keep away from Rain | | Batch Number | | Manufacturer |
| | Use by / Expiry date | | Temperature Limitation | | Read instructions for use | | Sterilised by Ethylene Oxide and single sterile barrier system |
| | Pyrogen free | | Keep away from sunlight | | Latex free product | | Does not contain DEHP phthalate |
| | Do not Resterilize | | For single use only | | Do not use if unit package is damaged | | Pressure |
| | The number of drops per milliliter | | Filter of liquid with pore size | | Medical Device | | Date and Country of Manufacture |
| | Volume of tubing | | Unique Device Identifier | | | | |

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