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## Infusion set IS-127 P® / Infusion set IS-127 P® 180 cm

CE 1014

<b>Catalogue No.:</b>	Medical device variants:
V606127-ND	Infusion set IS-127 P
V606128-ND	Infusion set IS-127 P 180 cm
<b>Product type:</b>	medical device
<b>Design:</b>	sterile, ethylene oxide sterilization
<b>Target patient:</b>	population of pediatric or adult patients
<b>Designated user:</b>	medical personnel with appropriate education, knowledge and practical experience in medical devices.
<b>Materials used:</b>	ABS, PE, PVC, PS, PP
<b>Compatibility:</b>	With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With plastic infusion bags, plastic or glass bottles. With infusion pharmaceuticals without any special requirements for filtration. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see <a href="http://www.sukl.cz">www.sukl.cz</a> or <a href="http://www.olecich.cz">www.olecich.cz</a>
<b>Incompatibility:</b>	With infusion pharmaceuticals incompatible with PVC (cytostatics) and with infusion pharmaceuticals sensitive to light.
<b>Contraindication:</b>	Do not use on any patient with a known allergic reaction to any of the product components.
<b>Adverse effects:</b>	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
<b>Warning:</b>	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Infusion bottle of contact with the infusion set until immediately before administration. During the operation, check the patient's condition and functionality sets. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
<b>Storage:</b>	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

### Intended use

The set is intended for therapeutic application of infusion solutions from bottles (bags) using gravitation or dosing infusion pumps:

INFUSOMAT P - made by B/Braun ARGUS A400, A404, A414, 707 - made by Argus

VOLUMED µVP5005 - made by Arcomed AG

The set shall only be used in connection with the type of pump specified!

### Instruction for use

- Under aseptic conditions, remove the set from the packing.
  - Close the hose regulation clamp.
  - Take off the cover of the piercing spike and the cover of the cone.
  - Inject the piercing spike upright into the plug of the bottle or bag. Hang the bottle or bag onto the stand.
  - Pressing the drop chamber fill it by 1/3, open the regulation clamp, open the plug of the air inlet opening and let the air out from the hose.
- Note: If using the set for transmission of solutions from a plastic bag, keep the air inlet opening closed.
- Close the regulation clamp.
  - Put the drop detector onto the drop chamber as described in the instructions of the manufacturer of the pump.
  - Open the hatch of the infusion pump and insert the stretched kit into the infusion pump.
  - Close the hatch of the infusion pump and open the regulation clamp.
  - Check if the solution does not flow through spontaneously.
  - Connect the set on the patient, set the required dosing rate in ml/h and switch on the infusion pump.

NOTE: The regulation clamp shall be completely open when the pump is in operation. After the dosing is finished, first close the regulation clamp and then take the kit out from the device.

### PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

### Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use
	Pyrogen free		Keep away from sunlight		Latex free product
	Do not Resterilize		For single use only		Do not use if unit package is damaged
	The number of drops per milliliter		Filter of liquid with pore size		Length of tubing
	Medical Device		Date of Manufacture		Volume of tubing

The instruction for use is also available in electronic form at [www.gama.cz](http://www.gama.cz). Free copies of other copies are available on request via e-mail: sales@gama.cz.

GAMA GROUP a.s., Betrieb Jimramov – 592 42 Jimramov  
 Sitz der Gesellschaft GAMA GROUP a.s., Mánesova 11/3b, České Budějovice 7, 370 01 České Budějovice, Tschechische Republik, Id.Nr.: 45797463, [www.gama.cz](http://www.gama.cz)

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VOLUMED µVP5005 - Hersteller Arcomed AG