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**EN Infusion set IS-127 P® M / Infusion set IS-127 P® 180 cm M**

**Catalogue No.:** V606127-M **Medical device variant:** Infusion set IS-127 P M **Catalogue No.:** V606128-M **Medical device variant:** Infusion set IS-127 P 180 cm M

**Functionality:** the manufacturer declares the functionality of the device when used in accordance with the instructions for use

**Duration of use:** „Short term“ - continuous use for the time maximum is 24 hours. Is dependent on the intended therapy according to the SPC of the drug/solution. Change the set according to national guidelines and/or hospital protocols.

**Design:** sterile, ethylene oxide sterilization

**Target patient:** Patients undergoing intravenous or intra-arterial infusion therapy.

**Designated user:** Physicians, anesthetists, nurses, and other healthcare professionals in healthcare settings who administer intravenous or intra-arterial infusion therapy.

**Materials used:** PVC, PE, ABS, PP, PA, membrane VERSAPORE RC, PS

**Compatibility:** With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With plastic infusion bags, plastic or glass bottles. With prescribed dosing pump. With infusion pharmaceuticals without any special requirements for filtration. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see [www.sukl.cz](http://www.sukl.cz) or [www.olecich.cz](http://www.olecich.cz)

**Incompatibility:** With infusion pharmaceuticals incompatible with PVC (cytostatics) and with infusion pharmaceuticals sensitive to light.

**Contraindication:** Do not use on any patient with a known allergic reaction to any of the product components.

**Adverse effects:** No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.

**Warning:** Do not use if package is open, damaged, guards are dropped or have expired. This may compromise the integrity or sterility of the device. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It can lead to contamination or impaired functionality, which can cause injury, illness, or death to the patient. Connect the infusion bottle to the infusion set just before administration. During the procedure, check the patient's condition and the functionality of the set. Use the set only in conjunction with the prescribed pump type!

**Storage and transport:** The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

**Intended use**

A single use medical device designed for therapeutic administration of infusion solutions from bottles (bags) by gravity or using infusion pumps: INFUSOMAT P – made by B/Braun a ARGUS A400, A404, A414, 707 – made by Argus.

**Instruction for use**

- Under aseptic conditions, remove the set from the packing.
- Close the hose regulation clamp.
- Take off the cover of the piercing spike and the cover of the cone.
- Inject the piercing spike upright into the plug of the bottle or bag. Hang the bottle or bag onto the stand.
- Pressing the drop chamber fill it by 1/3, open the regulation clamp, open the plug of the air inlet opening and let the air out from the hose.  
Note: If using the set for transmission of solutions from a plastic bag, keep the air inlet opening closed.
- Close the regulation clamp.
- Put the drop detector onto the drop chamber as described in the instructions of the manufacturer of the pump.
- Open the hatch of the infusion pump and insert the stretched kit into the infusion pump.
- Close the hatch of the infusion pump and open the regulation clamp.
- Check if the solution does not flow through spontaneously.
- Connect the set on the patient, set the required dosing rate in ml/h and switch on the infusion pump.

**NOTE:** The regulation clamp shall be completely open when the pump is in operation. After the dosing is finished, first close the regulation clamp and then take the kit out from the device.

**PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:**

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device.

Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

**Used symbols on packaging**

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide and single sterile barrier system
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Pressure
	The number of drops per milliliter		Filter of liquid with pore size		Length of tubing		Volume of tubing
	Medical Device		Date and Country of Manufacture		Unique Device Identifier		

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