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## Aspiration curette

<b>Catalogue No.:</b>	<b>Medical device variants:</b>
V646680-ND	Aspiration curette 6mm x 200mm
V646681-ND	Aspiration curette 7mm x 200mm
V646682-ND	Aspiration curette 8mm x 200mm
V646683-ND	Aspiration curette 9mm x 200mm
<b>Product type:</b>	medical device
<b>Design:</b>	sterile, ethylene oxide sterilization
<b>Target patient:</b>	population of girls and adult female patients
<b>Designated user:</b>	medical personnel with appropriate education, knowledge and practical experience in medical devices.
<b>Materials used:</b>	PVC
<b>Compatibility:</b>	With suction silicon or PVC tubes for secretion aspirators. The aspiration curette is not intended for drug delivery.
<b>Contraindication:</b>	Do not use on any patient with a known allergic reaction to any of the product components.
<b>Adverse effects:</b>	No adverse effects are known if the catheter is used in accordance with the instructions specified in the Warning and Use sections.
<b>Warning:</b>	Do not use, if the sterile packing is opened or damaged. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
<b>Storage:</b>	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

### Intended use

Curette is intended for gynaecological purposes - curettage.

### Instruction for use

1. Check the type labelling of the curette, usable life on the package and the package for damages.
2. Take the curette out of the cover.
3. Use the curette for the intended purpose.

### PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

### Used symbols on packaging

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Medical Device
	Date of Manufacture						

The instruction for use is also available in electronic form at [www.gama.cz](http://www.gama.cz). Free copies of other copies are available on request via e-mail: [sales@gama.cz](mailto:sales@gama.cz).

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