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Venous pressure meter

Catalogue No.: V646992-ND	Medical device variants: Venous pressure meter MŽT-01	Catalogue No.: V646992-01ND	Medical device variants: Venous pressure meter MŽT-01 MLL	Catalogue No.: V646993-ND	Medical device variants: Venous pressure meter MŽT-02
Product type: Design: Target patient: Designated user: Materials used: Compatibility:	medical device sterile, ethylene oxide sterilization population of pediatric or adult patients medical personnel with appropriate education, knowledge and practical experience in medical devices. PVC, PE, PC With infusion sets in conformity with standard ČSN EN ISO 8536-4 or ČSN EN ISO 8536-8. With intravenous catheters equipped with a connector with Luer thread in conformity with standard ČSN EN ISO 80369-7. The venous pressure meter is not intended for drug delivery.				
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components. Measurements under various conditions (change of position adjustment zero). After measurement the meter is disconnected.				
Adverse effects:	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.				
Warning:	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.				
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.				

Intended use


The venous pressure meter is used to determine the pressure in central veins of a patient. The meter is a disposable medical equipment to be used just once.

The venous pressure meter may be applied only when accessorized with an auxiliary device (not delivered as a part of the product) capable of providing:

- Vertical position of the related pressure-gage branch with a scale printed on it;
- Adjustment of scale zero point (0) so that the point corresponds with the height of jugular depression (fossa jugularis) in a lying patient.

Instruction for use

- Take the venous pressure meter out of the packing and place it onto the auxiliary device (see Figure 1).
- Remove the cover from the three-way valve and attach the branch tube to the infusion set. Take off the hose cover; the hose terminates in:
 - a positive cone (type: MŽT-01), or
 - a "Y" piece with a rubber stopper and a luer lock cone (type: MŽT-02).
- Turn the three-way valve to Position 2 (see Figure 3) and keep it there long enough for the pressure gage branch of the venous pressure meter to be filled up with the solution to the level of approx. 10 cm above the zero point of the scale (see Figure 1).
- Turn the three-way valve to Position 3 (see Figure 4) and read the solution level on the scale of the pressure gage branch (in cm). The solution level will be in proportion to the blood pressure in the patient's central veins.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use	Sterilised by Ethylene Oxide	
Pyrogen free		Keep away from sunlight	Latex free product	Does not contain DEHP phthalate	
Do not Resterilize		For single use only	Do not use if unit package is damaged	Medical Device	
Date of Manufacture					

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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