CS

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Double-ended piercing spike NO PVC

CE 1014

Catalogue No.:	V696439
Product type:	medical device
Design:	sterile, ethylene oxide sterilization
Target patient:	it is not intended for patients
Designated user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used:	ABS, PE
Compatibility:	The rubber stopper medication container (bottle). (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz
Warning:	Do not use, if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates it may lead to contamination or impairment functional.
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

The double-ended piercing spike is intended for the transfer, dilution and mixing of sterile solutions or drugs in a closed system between two containers (bottles), which are closed with a rubber stopper.

Instruction for use

- Under aseptic conditions remove the spike from the packaging.
- Remove the guards from the spike.
- Grasp the piercing spike's shaped holder and insert one end of the spike in one bottle's stopper and the other end in the other bottle's stopper.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Umbrella	Keep away from Rain	LOT	Batch Number	Factory	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book	Read instructions for use	Sterile EO	Sterilised by Ethylene Oxide
Pyrogen free	Pyrogen free	Sun	Keep away from sunlight	Latex	Latex free product	DEHP	Does not contain DEHP phthalate
Do not Resterilize	Do not Resterilize	Single use	For single use only	Do not use if unit package is damaged	Do not use if unit package is damaged	MD	Medical Device
Date of Manufacture	Date of Manufacture						

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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