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**EN Double-ended piercing spike NO PVC**

**Catalogue No.:** V696439  
**Functionality:** the manufacturer declares the functionality of the device when used in accordance with the instructions for use  
**Duration of use:** „Transient“ - continuous use for a maximum of 60 minute  
**Design:** sterile, ethylene oxide sterilization  
**Target patient:** it is not intended for patients  
**Designated user:** a healthcare worker with adequate education, knowledge and practical experience focused on this MD. The product is intended exclusively for professional use.  
**Materials used:** ABS, PE  
**Compatibility:** The rubber stopper medication container (bottle). With a closures for infusion bottles according to the ČSN EN ISO 8536-2 standard and with a plastic bag for intravenous infusions according to the ČSN EN ISO 15747 standard. With sterile solutions or lyophilized substances. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see [www.sukl.cz](http://www.sukl.cz) or [www.olecich.cz](http://www.olecich.cz)  
**Warning:** Do not use, if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions can lead to contamination or deterioration of the functionality of the MD.  
**Storage and transport:** The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

**Intended use**

Single-use medical device intended for the transfer, dilution and mixing of sterile solutions or drugs in a closed system between two containers (bottles), which are closed with a rubber stopper.

**Instruction for use**

- Under aseptic conditions remove the spike from the packaging.
- Remove the guards from the spike.
- Grasp the piercing spike's shaped holder and insert one end of the spike in one bottle's stopper and the other end in the other bottle's stopper.

**PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:**

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

**Used symbols on packaging**

REF	Catalogue Number	Umbrella	Keep away from Rain	LOT	Batch Number	Factory	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book with 'i'	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide and single sterile barrier system
Pyrogen	Pyrogen free	Sun with rays	Keep away from sunlight	Latex	Latex free product	PVC FREE	Does not contain PVC
Resterilize	Do not Resterilize	Single use	For single use only	Damaged	Do not use if unit package is damaged	MD	Medical Device
Date and Country	Date and Country of Manufacture	UDI	Unique Device Identifier				

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