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EN Transfer set SK


Catalogue No.: V606397-ND
Product type: medical device
Design: sterile, ethylene oxide sterilization
Target patient: population of pediatric or adult patients
Designated user: medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used: PVC
Compatibility: Syringe, single-use Luer-Slip or Luer-Lock. Funnel connector Charrière. The transfer set is not intended for drug delivery.
Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects: No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning: Do not use if the sterile packing is opened or damaged. Do not use after the expiry date.
 Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation.
 Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

The set is designed to transfer a flushing solution (flushing) where it is desirable to connect the syringe to a catheter or a flexible medical hose.

Instruction for use

1. Remove the transfer set out of the packing.
2. Connect the stepped connector to the catheter funnel end.
3. Attach the negative Female Luer Lock cone to the syringe.
4. Flush.

Procedure for safe disposal of medical devices after use:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Umbrella	Keep away from Rain	LOT	Batch Number	Factory	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide
Pyrogen	Pyrogen free	Sun	Keep away from sunlight	Latex	Latex free product	DEHP	Does not contain DEHP phthalate
Resterilize	Do not Resterilize	Single use	For single use only	Damaged	Do not use if unit package is damaged	MD	Medical Device
Date of Manufacture	Date of Manufacture						

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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