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EN Priming set of dialysis sets PS 003 M

Catalogue No.: V606003-M
Functionality: the manufacturer declares the functionality of the device when used in accordance with the instructions for use
Duration of use: „Short term“ - continuous use for the time maximum is 24 hours
Design: sterile, ethylene oxide sterilization
Target patient: A patient with renal failure.
Designated user: Doctors, nurses and other healthcare professionals in healthcare facilities who are qualified to perform hemodialysis.
Materials used: PS, PE, PVC, ABS, nylon
Compatibility: With dialysis sets equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With infusion plastic bags or plastic bottles. The Priming set of dialysis sets PS 003 is not intended for drug delivery.
Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.
Warning: Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage and transport: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

Intended use

A single use medical device intended for filling and rinsing dialysis sets and dialyzers before starting hemodialysis, or flushing blood from them when ending hemodialysis with physiological solution or, if necessary, topping up the patient with physiological solution.

The sets can be used at any type of dialyses units.

Instruction for use

- Under aseptic conditions, remove the set from the packing.
 - Close the pushers on both hoses.
 - Remove the cover of the piercer and the cone cover.
 - Stab the piercer upright into the stopper of the bottle or the bag. Hang the bottle or the bag on a stand.
 - Press the dropping chamber so that it fills up to 1/3, open the regulating pusher, open the stopper of the air-inlet opening and deaerate the hose.
- Note: If using the set to transfer the solution from a plastic bag, leave the air-inlet hole closed.
- Close the pusher.
 - Connect the arterial inlet of the dialysis unit to the short hose with the negative cone. Connect the long hose with the positive cone to the short hose in front of the pumping segment of the dialyses unit; leave the pusher closed.
 - Before starting the blood pump of the dialysis unit, open the short hose pusher.
 - After the dialyses unit and the dialyser are rinsed out, close the pusher.
 - The long hose remains connected to the dialysis unit during the whole dialysis; release the pusher if you want to fill up the physiological solution in course of the dialysis or in case of blood return from the dialysis set to the patient after the dialysis is finished.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use	Sterilised by Ethylene Oxide and single sterile barrier system	
Pyrogen free		Keep away from sunlight	Latex free product	Does not contain DEHP phthalate	
Do not Resterilize		For single use only	Do not use if unit package is damaged	Medical Device	
Date and Country of Manufacture		UDI	Unique Device Identifier		

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