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EN Priming set of dialysis sets PS 003

Catalogue No.: V606003-ND
Product type: medical device
Design: sterile, ethylene oxide sterilization
Target patient: it is not intended for patients
Designated user: medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used: PS, PA, PE, PVC
Compatibility: With dialysis sets equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With infusion plastic bags or plastic bottles. The Priming set of dialysis sets PS 003 is not intended for drug delivery.
Warning: Do not use, if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates it may lead to contamination or impairment of set functional.
Storage: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

The set is designated for rinsing of dialysis units using physiological solution from bottles (bags) or substitution solution from the dialysis apparatus.

The sets can be used at any type of dialyses units.

Instruction for use

1. Remove the set from the package.
 2. Close the pushers on both hoses.
 3. Remove the cover of the piercer and the cone cover.
 4. Stab the piercer upright into the stopper of the bottle or the bag. Hang the bottle or the bag on a stand.
 5. Press the dropping chamber so that it fills up to 1/3, open the regulating pusher, open the stopper of the air-inlet opening and deaerate the hose.
- Note: If using the set to transfer the solution from a plastic bag, leave the air-inlet hole closed.
6. Close the pusher.
 7. Connect the arterial inlet of the dialysis unit to the short hose with the negative cone. Connect the long hose with the positive cone to the short hose in front of the pumping segment of the dialyses unit; leave the pusher closed.
 8. Before starting the blood pump of the dialysis unit, open the short hose pusher.
 9. After the dialyses unit and the dialyser are rinsed out, close the pusher.
 10. The long hose remains connected to the dialysis unit during the whole dialysis; release the pusher if you want to fill up the physiological solution in course of the dialysis or in case of blood return from the dialysis set to the patient after the dialysis is finished.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide
Pyrogen free		Keep away from sunlight	Latex free product	Does not contain DEHP phthalate	
Do not Resterilize		For single use only	Do not use if unit package is damaged	MD	Medical Device
Date of Manufacture					

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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