## CS Prodlužovací hadička PVC, 150 cm M

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## extension PVC tube 150 cm M

**Catalogue No.:** V606502-M  
**Functionality:** the manufacturer declares the functionality of the device when used in accordance with the instructions for use  
**Duration of use:** „Short term“. If blood, blood derivatives, fat emulsions are administered or if the line had to be disconnected, continuous use is for a maximum of 24 hours, in case of administration of the same type of drug, use is for a maximum of 72 hours, in a closed system terminated with a needle-free inlet the interval can be extended up to 7 days. Immediate replacement is necessary in case of contamination of the line and in case of incompatibility. Replacement according to national guidelines and/or procedures of the healthcare facility. Depends on the intended treatment according to the summary of product characteristics of the drug/solution.  
**Design:** sterile, ethylene oxide sterilization  
**Target patient:** Patients undergoing intravenous or intra-arterial infusion therapy where an extension or flexible connection between the infusion system and the intravenous/intra-arterial catheter is required.  
**Designated user:** Doctors, anesthesiologists, nurses and other healthcare workers in healthcare facilities who perform intravenous or intra-arterial infusion therapy.  
**Materials used:** PVC, PE, ABS  
**Compatibility:** With medical devices equipped with standard Luer connectors in conformity with standards ČSN EN ISO 80369-7. Extension PVC tubes are designed to transfer liquid or i.v. administration. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see [www.sukl.cz](http://www.sukl.cz) or [www.olecich.cz](http://www.olecich.cz).  
**Contraindication:** Do not use on any patient with a known allergic reaction to any of the product components.  
**Adverse effects:** No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.  
**Warning:** Use up to a pressure of max. 2 bar. Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.  
**Storage and transport:** The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

### Intended use

Single use medical devices intended as connecting tubes for linear dispensers for the administration intravenous or intra-arterial infusion therapy. Tubing terminations on both ends of the Luer Lock cone ensure a secure connection to a standard syringe designed for a linear dispenser and to an intravenous or intra-arterial catheter. Furthermore, this tubing can be used as a standard extension tubing to connect an infusion set to an intravenous needle. The connection connector of the extension tube are compatible with Luer Lock connectors, which meet the requirements of the ČSN EN ISO 80369-7 standard.

### Instruction for use

- Under aseptic conditions remove the tube from the packaging.
- Remove the guards from the connecting beams.
- Attach the tube's negative cone to the positive cone of the set (on the syringe) and the positive cone to the catheter.

#### PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

### Used symbols on packaging

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide and single sterile barrier system
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Medical Device
	Date and Country of Manufacture		Unique Device Identifier		Pressure		

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