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EN Flow regulator infusion

Catalogue No.:	V606145-ND
Product type:	medical device
Design:	sterile, ethylene oxide sterilization
Target patient:	population of pediatric or adult patients
Designated user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used:	PVC, PE, PP, ABS, PC, POM
Compatibility:	With gravity infusion set, which is in conformity with standard ČSN EN ISO 8536-4. With intravenous cannula with Luer Lock connector. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects:	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning:	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

Disposable medical device for gravity infusion flow control. The flow regulator infusion may be used for saline of viscosity from 0.9% NaCl up to 40% of glucose solution. It is universal for connection to all types of the infusion sets. It is recommended for use with 20G IV cannulas. Adjustable infusion flow rate from 10 to 300 ml/hour.

Instruction for use

- Under aseptic conditions, unpack the flow regulator infusion and remove protective covers from the connection cones.
 - Connect the negative connection cone to the infusion gravity set.
 - Inject the vented spike of the infusion set perpendicularly to the bottle or bag cap and close the regular clamp.
 - Hook the bottle or bag on the IV stand so that the height between the solution level in the drip chamber and output (end) of the I.V. catheter is 80cm.
 - Press the drip chamber to fill in the same with 1/3 of the solution, open the regular clamp, open plug of the aeration hole and purge air from the hose.
 - Close the regular clamp and rotate the regulator (from OPEN to OFF position) to set the required flow rate in ml/hour.
 - Open the regular clamp pusher and check configured flow rate by counting the drops. If needed, adjust height of the bottle or bag to increase or decrease the flow.
 Example of the drop counting:
 200 ml/h with drip chamber for 20 drops = 1 ± 0.1 ml means 66 drops per minute
 100 ml/h with drip chamber for 20 drops = 1 ± 0.1 ml means 33 drops per minute
 - Close the regular clamp and connect the positive connecting cone to the IV cannula of the patient.
 - Open the regular clamp and start feeding of the infusion.
- ⚠ Use for gravity infusions only. Do not use for feeding of blood, blood derivatives or solutions of viscosity over 40%. Check the configured flow rate by counting the drops during the infusion. To change the flow rate, adjust infusion solution bottle or bag height. Height between the solution level in the drip chamber and output (end) of the I.V. catheter should be 80cm. Change to height of the infusion solution by 10cm increases flow rate by about 13%. Use only under flow rate supervision, particularly for the infusions where too high or too low flow rate may be hazardous. The flow regulator is delivered always in open position.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:
 Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number		Keep away from Rain	LOT	Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Attention
MD	Medical Device		Date of Manufacture				

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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