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**EN Flow regulator infusion**

**Catalogue No.:** V606145-ND  
**Functionality:** the manufacturer declares the functionality of the device when used in accordance with the instructions for use  
**Duration of use:** „Short term“ - continuous use for the time maximum is 24 hours. Is dependent on the intended therapy according to the SPC of the drug/solution. Change the regulator according to national guidelines and/or hospital protocols.  
**Design:** sterile, ethylene oxide sterilization  
**Target patient:** population of pediatric or adult patients  
**Designated user:** a healthcare worker with adequate education, knowledge and practical experience focused on this MD. The product is intended exclusively for professional use.  
**Materials used:** PE, PVC, ABS, POLYISOPRENE, PP, PC, POM  
**Compatibility:** With gravity infusion set, which is in conformity with standard ČSN EN ISO 85363-4. With an intravenous cannula with a Luer Lock connector that complies with the ISO 80369-7 standard. With infusion drugs that do not have special filtration requirements intended for intravenous administration. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see [www.sukl.cz](http://www.sukl.cz) or [www.olecich.cz](http://www.olecich.cz)  
**Contraindication:** Do not use on any patient with a known allergic reaction to any of the product components.  
**Adverse effects:** No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.  
**Warning:** Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date.  
 Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.  
**Storage and transport:** The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation.  
 Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action.  
 The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them.  
 The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

**Intended use**

Disposable medical device for gravity infusion flow control. The flow regulator infusion may be used for saline of viscosity from 0.9% NaCl up to 40% of glucose solution. It is universal for connection to all types of the infusion sets. It is recommended for use with 20G IV cannulas. Adjustable infusion flow rate from 10 to 350 ml/hour.

**Instruction for use**

- Under aseptic conditions, unpack the flow regulator infusion and remove protective covers from the connection cones.
  - Connect the negative connection cone to the infusion gravity set.
  - Inject the vented spike of the infusion set perpendicularly to the bottle or bag cap and close the regular clamp.
  - Hook the bottle or bag on the IV stand so that the height between the solution level in the drip chamber and output (end) of the I.V. catheter is 1 m.
  - Press the drip chamber to fill in the same with 1/3 of the solution, open the regular clamp, open plug of the aeration hole and purge air from the hose.
  - Close the regular clamp and rotate the regulator (from OPEN to OFF position) to set the required flow rate in ml/hour.
  - Open the regular clamp pusher and check configured flow rate by counting the drops. If needed, adjust height of the bottle or bag to increase or decrease the flow.  
 Example of the drop counting:  
 200 ml/h with drip chamber for 20 drops = 1 ± 0.1 ml means 66 drops per minute  
 100 ml/h with drip chamber for 20 drops = 1 ± 0.1 ml means 33 drops per minute
  - Close the regular clamp and connect the positive connecting cone to the IV cannula of the patient.
  - Open the regular clamp and start feeding of the infusion.
- ⚠ Use for gravity infusions only. Do not use for feeding of blood, blood derivatives or solutions of viscosity over 40%. Check the configured flow rate by counting the drops during the infusion. To change the flow rate, adjust infusion solution bottle or bag height. Height between the solution level in the drip chamber and output (end) of the I.V. catheter should be 1 m. Change to height of the infusion solution by 10cm increases flow rate by about 13%. Use only under flow rate supervision, particularly for the infusions where too high or too low flow rate may be hazardous. The flow regulator is delivered always in open position.

**Dosing accuracy:** at flow rates up to 20 ml / h, dosing accuracy is within a tolerance range of ± 35%, at flow rates of 30 ml / h ± 25%, for flow rates higher than 30 ml / h ± 15%.

**PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:**  
 Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

**Used symbols on packaging**

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide and single sterile barrier system
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Attention
	Medical Device		Date and Country of Manufacture		Unique Device Identifier		

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