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
**FUNCTIONALITY:** the manufacturer declares the functionality of the device when used in accordance with the instructions for use.

**DURATION OF USE:** „Transient“ - continuous use for a maximum of 60 minute.

**DESIGN:** sterile, ethylene oxide sterilization.

**TARGET PATIENT:** a patient with advanced stage cancer that does not respond to standard systemic treatment

**DESIGNATED USER:** a qualified surgeon who is experienced in oncological surgery and certified in peritoneal surgery

**MATERIALS USED:** PVC, PE, ABS, POM, PP, PTFE, PE film, Membrane Versapore, EVA, SBC

**PURPOSE OF USE:** The set is single-use medical device intended for hyperthermic lavage of the abdominal cavity using a cytostatic solution. The set is designed for use with the device manufactured by SKALA-Medica s.r.o (SKALA LA-1).

**COMPATIBILITY:** The set is compatible with the SKALA LA-1 device (manufacturer SKALA-Medica s.r.o.) and with a plastic infusion bag. The set is intended for the administration of solutions and cytostatics used in oncology for intravenous administration. The set is intended for the administration of solutions and cytostatics used in oncology for intravenous administration. Check the (in)compatibility of the drug with PVC and ABS materials in the summary of product characteristics (SPC), which is listed on the SÚKL information portal in the drug database on the website www.sukl.cz or www.olecich.cz..

**DESCRIPTION**

SKALA OK3 Perfusion Set comprises 7 parts: SKALA OK3/1 Main tube set, SKALA OK3/2 Exchanger tubes, SKALA OK3/3 Temperature measuring tube set part, SKALA OK3/4 Patient parts, SKALA OK3/5 Additional filter, SKALA OK3/6 Container 3000 ml, SKALA OK3/7 Protective surgical sleeve. V606501-M has all the parts. V606501-01M, V606501-02M and V606501-03M does not have Protective surgical sleeve. V606501-03 does not have Additional filter.

Depending on the method of the lavage performed, it is possible to choose a set type with an adjusted patient part for the specific method:

The SKALA OK3 perfusion set incorporates a patient part intended to perform open lavage.

The SKALA OK3-V1 perfusion set variant of the lavage set incorporates a patient part intended for performing closed lavage

The SKALA OK3-V2 perfusion set variant of the lavage set incorporates a patient part intended for the treatment of ascites; all other parts constituting the set remain the same

**NOTICE**

Do not use the set in case of non-conformity with the type marking on the package, damage to the package or dropped guards. Check the package for completeness before use.

**STORAGE AND TRANSPORT**

The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation.

Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action.

The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

**INSTRUCTIONS FOR USE**


- Remove the bag “SKALA OK3/6 – Container” from the package and hang on a hook intended for this purpose.
- Remove the part “SKALA OK3/2 - Exchanger tubes” from the package.
- Insert the part “SKALA OK3 - Exchanger tubes” to the exchanger container. Leave the connectors (1) and (2) hanging freely over the edge of the container and cover the container.
- Remove the part “SKALA OK3/1 – Main tube set” from the package.
- Insert the drip chamber (5) in the holder on the machine marked HOLDER.
- Insert the double segment of the pump (17) in the path of the peristaltic pump so that the inlet part of the segment IN and the outlet part OUT match the marks on the machine.
- Insert a part of the inlet tube from the patient (15) under the automatic valve PAT.OUT and in the path of the blood leak detector. Close the pinch-cock (14) on this tube.
- Insert the piercing tool (10) in the container outlet. Insert a part of the adjacent tube (12) under the automatic valve CONT.OUT.
- Connect the connector of the separating filter (blue cover) (16) to the inlet pressure IN meter in the back of the machine.
- Connect the connector of the separating filter (red cover) (7) to the outlet pressure OUT meter in the back of the machine. Close the pinch-cock (6) on the spare separating filter.
- Connect a syringe to the separating filter (3) and close the pinch-cock (4).
- Close the pinch-cocks (8).
- Connect the connectors (9) to the inlet tubes of the exchanger (2).
- Remove the part “SKALA OK3/3 – Temperature measuring tube set part” from the package.
- Connect the connectors (24) to the outlet tubes of the exchanger (1).
- Insert a part of the outlet tube leading to the patient (21) under the automatic valve PAT.IN and close the adjacent pinch-cock (20).
- Connect the connector (18) to the container inlet. Insert a part of the adjacent tube (22) under the automatic valve CONT.IN.
- \*) \*\*) Remove the part “SKALA OK3/7 – Protective surgical sleeve” from the package and use it to cover the operation wound.
- Remove the part “SKALA OK3/4 – Patient parts” from the package.
- Insert the double ends of both parts in the patient’s abdominal cavity.
- Connect the patient outlet connector (25) to the connector (13).
- Connect the patient inlet connector (28) to the connector (19).
- Before filling the patient, open the pinch-cocks (20) and (14). \*\*\*)
- After use, products should be disposed of in line with prescribed procedures for disposal of medical (infectious) waste in accordance with local regulations.

\*) Steps 18. – 20. shall be performed by operating personnel.

\*\*) Only in case of open lavage.

\*\*\*) If you use option sets SKALA OK3-V1 and SKALA OK3-V2 enclose the filling patient clamp (31), remove the protective cover from the connector (29) and (32) and after soaking the patient to close the clamp for venting (30). Before the extraction solution from the patient to open the clamp for venting.

**NOTE** – if the solution passes through the separating filter at the connector (7), close the adjacent pinch-cock, disconnect the connector (7) from the pressure meter and instead connect the reserve branch (R1) to the pressure meter.  
 After that, release the pinch-cock (6).  
 – in case the filter is clogged, connect the part SKALA OK3/5 Additional filter to the connectors (R2) and release the pinch-locks (8).

**SKALA OK3 PERFUSIONSET M / SKALA OK3-V1 PERFUSIONSET M / SKALA OK3-V2 PERFUSIONSET M / SKALA OK3-V3 PERFUSIONSET M**

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