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

	= Sterilised using ethylene oxide		= Do Not Reuse		= Batch Code
	= Non-pyrogenic		= Attention, see instructions for use		= Use by
	= Catalogue Number		= Manufacturer		= contains di (2-ethylhexyl)-phthalate (DEHP)
	= Do not Resterilize		= Date of Manufacture		

PURPOSE OF USE:
 The set is intended for hyperthermic lavage of the abdominal cavity using a cytotatic solution. This set is designed to be used with the SKALALA-1 machine (manufacturer SKALA-Medica s.r.o.).


COMPATIBILITY:
 The set is intended for the administration of solutions and cytotatics used in oncology for intravenous administration. (In) compatibility of the PVC and ABS drug with the Summary of Product Characteristics (SPC) given on the information portal SÚKL in the drug database at www.sukl.cz or www.olecich.cz. The set is compatible with the SKALA LA-1 (manufacturer SKALA-Medica s.r.o.) and a plastic infusion bag.

DESCRIPTION:
 SKALA OK3 Perfusion Set comprises 7 parts: SKALA OK3/1 Main tube set, SKALA OK3/2 Exchanger tubes, SKALA OK3/3 Temperature measuring tube set part, SKALA OK3/4 Patient parts, SKALA OK3/5 Additional filter, SKALA OK3/6 Container 3000 ml, SKALA OK3/7 Protective surgical sleeve.
 Depending on the method of the lavage performed, it is possible to choose a set type with an adjusted patient part for the specific method:
 The SKALA OK3 perfusion set incorporates a patient part intended to perform open lavage.
 The SKALA OK3-V1 variant of the lavage set incorporates a patient part intended for performing closed lavage
 The SKALA OK3-V2 variant of the lavage set incorporates a patient part intended for the treatment of ascites; all other parts constituting the set remain the same


NOTICE:
 Do not use the set in case of non-conformity with the type marking on the package, damage to the package or dropped guards. Check the package for completeness before use.

STORAGE:
 The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

- INSTRUCTIONS FOR USE:**
- Remove the bag "SKALA OK3/6 – Container" from the package and hang on a hook intended for this purpose.
 - Remove the part "SKALA OK3/2 - Exchanger tubes" from the package.
 - Insert the part "SKALA OK3 - Exchanger tubes" to the exchanger container. Leave the connectors (1) and (2) hanging freely over the edge of the container and cover the container.
 - Remove the part "SKALA OK3/1 – Main tube set" from the package.
 - Insert the drip chamber (5) in the holder on the machine marked HOLDER.
 - Insert the double segment of the pump (17) in the path of the peristaltic pump so that the inlet part of the segment IN and the outlet part OUT match the marks on the machine.
 - Insert a part of the inlet tube from the patient (15) under the automatic valve PAT.OUT and in the path of the blood leak detector. Close the pinch-cock (14) on this tube.
 - Insert the piercing tool (10) in the container outlet. Insert a part of the adjacent tube (12) under the automatic valve CONT.OUT.
 - Connect the connector of the separating filter (blue cover) (16) to the inlet pressure IN meter in the back of the machine.
 - Connect the connector of the separating filter (red cover) (7) to the outlet pressure OUT meter in the back of the machine. Close the pinch-cock (6) on the spare separating filter.
 - Connect a syringe to the separating filter (3) and close the pinch-cock (4).
 - Close the pinch-cocks (8).
 - Connect the connectors (9) to the inlet tubes of the exchanger (2).
 - Remove the part "SKALA OK3/3 – Temperature measuring tube set part" from the package.
 - Connect the connectors (24) to the outlet tubes of the exchanger (1).
 - Insert a part of the outlet tube leading to the patient (21) under the automatic valve PAT.IN and close the adjacent pinch-cock (20).
 - Connect the connector (18) to the container inlet. Insert a part of the adjacent tube (22) under the automatic valve CONT.IN.
 - *, **) Remove the part "SKALA OK3/7 – Protective surgical sleeve" from the package and use it to cover the operation wound.
 - Remove the part "SKALA OK3/4 – Patient parts" from the package.
 - Insert the double ends of both parts in the patient's abdominal cavity.
 - Connect the patient outlet connector (25) to the connector (13).
 - Connect the patient inlet connector (28) to the connector (19).
 - Before filling the patient, open the pinch-cocks (20) and (14). ***)
 - After use, products should be disposed of in line with prescribed procedures for disposal of medical (infectious) waste in accordance with local regulations.


SKALA OK3/1 Main tube set

SKALA OK3/2 Exchanger tubings

*) Steps 18. - 20. shall be performed by operating personnel.
 **) Only in case of open lavage.
 ***) If you use option sets SKALA OK3-V1 and SKALA OK3-V2 enclose the filling patient clamp (31), remove the protective cover from the connector (29) and (32) and after soaking the patient to close the clamp for venting (30). Before the extraction solution from the patient to open the clamp for venting.

NOTE – if the solution passes through the separating filter at the connector (7), close the adjacent pinch-cock, disconnect the connector (7) from the pressure meter and instead connect the reserve branch (R1) to the pressure meter.
 After that, release the pinch-cock (6).
 – in case the filter is clogged, connect the part SKALA OK3/5 Additional filter to the connectors (R2) and release the pinch-locks (8).

SKALA OK3/3 Temperature measuring tube set part
 SKALA OK3/4 Patient's tube set segment for open lavage
 SKALA OK3/4-V1 Patient's tube set segment for closed lavage
 SKALA OK3/4-V2 Patient's tube set segment for treating ascites

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