CS **Spojovací hadička M**

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EN **Connecting tubing M**

Catalogue No.: V606402-M
V606404-M

Medical device variants:
Connecting Tubing 2-way M
Connecting Tubing 4-way M

Functionality: the manufacturer declares the functionality of the device when used in accordance with the instructions for use

Duration of use: „Short term“- continuous use for the time maximum is 24 hours

Design: sterile, ethylene oxide sterilization

Target patient: Patients undergoing intravenous or intra-arterial infusion therapy where multiple infusion sets need to be connected to an intravenous/intra-arterial catheter.

Designated user: Doctors, anesthesiologists, nurses and other healthcare workers in healthcare facilities who perform intravenous or intra-arterial infusion therapy.

Materials used: PVC, PE, PS, PP

Compatibility: With medical devices equipped with standard Luer connectors in conformity with standards ČSN EN ISO 80369-7. Connecting Tubing are designed to transfer liquid or i.v. administration or transfer of blood. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz.

Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.

Adverse effects: No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.

Warning: Use up to a pressure of max. 2 bar. Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.

Storage and transport: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

Intended use

Single use medical devices designed to connect multiple infusion sets with other components such as i.v. cannula. One or more infusion solutions or drug solutions can be administered simultaneously through the inlets of the connecting tubes, which form the Luer Lock connection cones. The connecting tubes are designed for pressure and gravity systems up to 2 bar pressure and are compatible with Luer Lock cones that meet the requirements of the ČSN EN ISO 80369-7 standard.

Instruction for use

- Under aseptic conditions remove the tube from the packaging.
- Remove the caps from the connecting beams.
- Attach the tube's negative cone to the infusion set (or syringe) and the positive cone to the intravenous cannula.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer	
	Use by / Expiry date	Temperature Limitation		Read instructions for use	Sterilised by Ethylene Oxide and single sterile barrier system	
	Pyrogen free	Keep away from sunlight		Latex free product	Does not contain DEHP phthalate	
	Do not Resterilize			Do not use if unit package is damaged	Pressure	
MD	Medical Device		VOL	Volume of tubing	UDI	Unique Device Identifier

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