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**EN Connection tubing with needleless port M**

**Catalogue No.:** V686424-M  
**Functionality:** the manufacturer declares the functionality of the device when used in accordance with the instructions for use  
**Duration of use:** „Short term“ - continuous use for the life of the IV cannula to which it is connected. However, the tubing should be replaced after a maximum of 7 days or 100 connections.  
**Design:** sterile, ethylene oxide sterilization  
**Target patient:** Patients undergoing intravenous or intra-arterial infusion therapy where an extension or flexible connection between the infusion system and the intravenous/intra-arterial catheter is required.  
**Designated user:** Doctors, anesthesiologists, nurses and other healthcare workers in healthcare facilities who perform intravenous or intra-arterial infusion therapy.  
**Materials used:** PC, PE, PVC, ABS, silicon  
**Compatibility:** With medical devices equipped with standard Luer connectors in conformity with standards ČSN EN ISO 80369-7. Connection tubing with needleless port are designed to transfer liquid or i.v. administration or transfer of blood. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see [www.sukl.cz](http://www.sukl.cz) or [www.olecich.cz](http://www.olecich.cz).  
**Contraindication:** Do not use on any patient with a known allergic reaction to any of the product components.  
**Adverse effects:** No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.  
**Warning:** Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date.  
 Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.  
**Storage and transport:** The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation.  
 Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action.  
 The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

**Intended use**

The connection tubing with needle-free entry is a disposable medical device designed to connect i.v. cannulas with other infusion systems with Luer ends (syringes, infusion sets, etc.). It can also be used to extend or connect any infusion set that has a Luer end. The connecting tube with needle-free entry ensures a closed system and thus microbiological safety for the period of use of the i.v. of the cannula to which it is connected, but its replacement must always be done after a maximum of 7 days or after 100 connections. The needle-free entrance with which the tube ends allows multiple and safe access (without the use of a needle) for direct injection or infusion. The connecting tube with needle-free entry is intended for pressure and gravity systems up to 2 bar pressure.

**Instruction for use**

- Under aseptic conditions remove the tube from the packaging.
- Remove the cover from the connecting cone.
- Before application, wipe or spray the surface of the silicone diaphragm (of the needleless port) with an alcoholic disinfectant.
- Wait for no less than 30 seconds until the disinfectant solution dries out.
- Attach the positive cone of the tubing (Male Luer Lock) by a clockwise twisting motion to an intravenous cannula so as to ensure a safe connection.
- Attach the negative cone of the needleless port by a clockwise twisting motion to an infusion system with a Luer ending so as to ensure a safe connection.
- If necessary, close the flow tube using clamp.
- Detach the tubing by a counter clockwise twisting motion.

**PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:**

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

**Used symbols on packaging**

REF	Catalogue Number	Umbrella	Keep away from Rain	LOT	Batch Number	Factory	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide and single sterile barrier system
X	Pyrogen free	Sun	Keep away from sunlight	Latex	Latex free product	PHT/DEHP	Does not contain DEHP phthalate
2	Do not Resterilize	X	For single use only	X	Do not use if unit package is damaged	P	Pressure
MD	Medical Device	Calendar	Date and Country of Manufacture	VOL	Objem hadičky	UDI	Unique Device Identifier

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