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**EN Connection tubing with needleless port**

**Catalogue No.:** V686424-ND  
**Product type:** medical device  
**Design:** sterile, ethylene oxide sterilization  
**Target patient:** population of pediatric or adult patients  
**Designated user:** medical personnel with appropriate education, knowledge and practical experience in medical devices.  
**Materials used:** PC, PE, PVC, ABS  
**Compatibility:** With medical devices equipped with standard Luer connectors in conformity with standards ČSN EN ISO 80369-7. Connection tubing with needleless port are designed to transfer liquid or i.v. administration or transfer of blood. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see [www.sukl.cz](http://www.sukl.cz) or [www.olecich.cz](http://www.olecich.cz).  
**Contraindication:** Do not use on any patient with a known allergic reaction to any of the product components.  
**Adverse effects:** No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.  
**Warning:** Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.  
**Storage:** The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

**Intended use**

Connection tubing with needleless port is designed to connect a peripheral or central intravenous catheter to other infusion systems with Luer ending (syringes, infusion sets, etc.). It can also be used to extend or connect any infusion system with Luer ending. The maximum time of serviceability depends on the durability of the catheter to which the tubing is connected, however it has to be replaced after 7 days of use or after 100 connections at most. The advantage of this product is the closed system that provides microbiological safety for the period of 7 days and the possibility of multiple and safe application (up to 100 applications). Connection tubing with needleless port is designed for both pressure and gravity systems up to the pressure of 2 bar. Connecting cones of the connection tubing are compatible with the Luer cones.

**Instruction for use**

- Under aseptic conditions remove the tube from the packaging.
- Remove the cover from the connecting cone.
- Before application, wipe or spray the surface of the silicone diaphragm (of the needleless port) with an alcoholic disinfectant.
- Wait for no less than 30 seconds until the disinfectant solution dries out.
- Attach the positive cone of the tubing (Male Luer Lock) by a clockwise twisting motion to an intravenous cannula so as to ensure a safe connection.
- Attach the negative cone of the needleless port by a clockwise twisting motion to an infusion system with a Luer ending so as to ensure a safe connection.
- If necessary, close the flow tube using clamp.
- Detach the tubing by a counter clockwise twisting motion.

**PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:**

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

**Used symbols on packaging**

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use		Sterilised by Ethylene Oxide
Pyrogen free		Keep away from sunlight	Latex free product		Does not contain DEHP phthalate
Do not Resterilize		For single use only	Do not use if unit package is damaged		Pressure
MD	Medical Device	Date of Manufacture	VOL	Objem hadičky	

The instruction for use is also available in electronic form at [www.gama.cz](http://www.gama.cz). Free copies of other copies are available on request via e-mail: [sales@gama.cz](mailto:sales@gama.cz).

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