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## Transfusion set for blood bags TSV-221

CE 1014

Catalogue No.:	V606221-ND
Product type:	medical device
Design:	sterile, ethylene oxide sterilization
Target patient:	population of pediatric or adult patients
Designated user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used:	PE, PVC, ABS, PS, PP
Compatibility:	With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With blood bags. With blood bags. The transfusion set for blood bags TSV-221 is not intended for drug delivery.
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects:	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning:	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

### Intended use

The set is designed to transfer full blood from a plastic bag to the patient's vein.

### Instruction for use

- Check the type designation on the set packing, inspect the packing if undamaged and make sure no cover of the set dropped off. Do not use the set if there is a discrepancy between the set and its type designation or in case the packing is damaged, or when the protective covers dropped off.
- Under aseptic conditions, remove the set from the packing. Close the control clamp (by moving the wheel to the narrow end); remove the cap from the plastic needle and drive the needle to the appropriate outlet of the blood bag.
- Hang the bag on a rack and press the drip chamber repeatedly so that it fills up with blood to approx. one half.
- Remove the cover from the connection cone and fill the hose full of blood by partly opening the control clamp. Then close the clamp.
- Now the set is ready to perform the transfusion. Start the transfusion in conformity to the relevant in-house regulation (as established in the hospital).

### PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

### Used symbols on packaging

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Filter of liquid with pore size
	The number of drops per milliliter		Medical Device		Date of Manufacture		

The instruction for use is also available in electronic form at [www.gama.cz](http://www.gama.cz). Free copies of other copies are available on request via e-mail: sales@gama.cz.

 GAMA GROUP a.s., Betrieb Jimramov – 592 42 Jimramov  
Sitz der Gesellschaft GAMA GROUP a.s., Mánesova 11/3b, České Budějovice 7, 370 01 České Budějovice, Tschechische Republik, Id.Nr.: 45797463, [www.gama.cz](http://www.gama.cz)

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