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Transfusion set for blood bags TSV-222

CE 1014

Catalogue No.:	V606222-ND
Product type:	medical device
Design:	sterile, ethylene oxide sterilization
Target patient:	population of pediatric or adult patients
Designated user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used:	PE, ABS, PVC, PS, PP, STEEL
Compatibility:	With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With blood bags. The transfusion set for blood bags TSV-222 is not intended for drug delivery.
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects:	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning:	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

The set is designed to dilute packed erythrocytes contained in a plastic bag and convey them, when diluted, to the patient's vein. The set can also be used to perform a simple plasmapheresis.

Instruction for use

- Check the type designation on the set packing, inspect the packing if undamaged and make sure no cover of the set dropped off. Do not use the set if there is a discrepancy between the set and its type designation or a in case the packing is damaged, or when the protective covers dropped off.
- The first item to be taken out of the packing is the aerator under aseptic conditions; remove the needle cover and drive the needle through stopper into a bottle containing saline.
- Then remove the set proper from the packing; close the control clamp and press closed also the two pinchcocks on the hoses upstream (above) of the drip chamber.
- Take off the cover of the metal needle (the one with the lesser diameter) and drive the needle through stopper into a bottle containing saline. Hang the bottle on a rack, taking care not to wet the aerator filter. With the bottle hung, fix the aerator hose firm so that the filter is above the saline level.
- Take off the cover of the plastic needle and drive it to the outlet of the relevant plastic bag containing the packed erythrocytes.
- Open both pinchcocks and transfer a sufficient volume of the saline from the bottle to the bag; close both the pinchcocks. Move the bag so that its content becomes thoroughly mixed; then hang the bag on a rack.
- Open the pinchcock at the outlet of the bag and pressed the drip chamber repeatedly so that it fills up with blood to approx. one half.
- Remove the cover of the connection cone and fill the hose full of blood by partly opening the control clamp. Then close the clamp.
- Now the set is ready to perform the transfusion. Start the transfusion in conformity to the relevant in-house regulation (established in the hospital).

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
	Use by / Expiry date			Read instructions for use	STERILE EO Sterilised by Ethylene Oxide
	Pyrogen free			Latex free product	Does not contain DEHP phthalate
	Do not Resterilize			Do not use if unit package is damaged	Filter of liquid with pore size 15 µm
	The number of drops per milliliter			Date of Manufacture	

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

GAMA GROUP a.s., Betrieb Jimramov – 592 42 Jimramov
 Sitz der Gesellschaft GAMA GROUP a.s., Mánesova 11/3b, České Budějovice 7, 370 01 České Budějovice, Tschechische Republik, Id.Nr.: 45797463, www.gama.cz

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