CS

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EN

Transfusion set TS-203 M / Transfusion set TS-204 M

Catalogue No.: V606203-M
 V606204-M

Medical device variants:
 Transfusion set TS-203 M
 Transfusion set TS-204 M

Functionality: the manufacturer declares the functionality of the device when used in accordance with the instructions for use

Duration of use: „Short term“ - continuous use for the time maximum is 24 hours

Design: sterile, ethylene oxide sterilization

Target patient: Patients undergoing intravenous or intra-arterial transfusion therapy.

Designated user: Doctors, anesthesiologists, nurses, and other healthcare workers in healthcare facilities who perform transfusion therapy.

Materials used: PE, ABS, PVC, PS, PP, PA, membrane VERSAPORE RC

Compatibility: With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With blood bags.

Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.

Adverse effects: No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.

Warning: Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.


Storage and transport: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

Intended use

Single-use medical device intended for therapeutic transferring of preserved blood, basic blood components, blood substitutes etc. from bottles or plastic bags into the patient.

Instruction for use

- Under aseptic conditions, remove the set from the packing.
- Close the hose regulation clamp.
- Take off the cover of the piercing spike and the cover of the cone.
- Inject the piercing spike upright into the plug of the transfusion bag. Hang the transfusion bag onto the stand.
- Pressing the drop chamber fill it by 1/3, open the regulation clamp, open the plug of the air inlet opening and let the air out from the hose.
- Close the regulation clamp.
- Connect the set on the patient, open the regulation clamp and start the transfusion.

Procedure for safe disposal of medical devices after use:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

| | | | | | | | |
|--|------------------------------------|--|-------------------------|--|---------------------------------------|--|--|
| | Catalogue Number | | Keep away from Rain | | Batch Number | | Manufacturer |
| | Use by / Expiry date | | Temperature Limitation | | Read instructions for use | | Sterilised by Ethylene Oxide and single sterile barrier system |
| | Pyrogen free | | Keep away from sunlight | | Latex free product | | Does not contain DEHP phthalate |
| | Do not Resterilize | | For single use only | | Do not use if unit package is damaged | | Filter of liquid with pore size TS-204 |
| | The number of drops per milliliter | | Medical Device | | Date and Country of Manufacture | | Filter of liquid with pore size TS-201 and TS-203 |
| | Unique Device Identifier | | | | | | |

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