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EN Transfusion set TS-201 / Transfusion set TS-203 / Transfusion set TS-204


Catalogue No.: V606201-ND, V606203-ND, V606204-ND
Medical device variants: Transfusion set TS-201, Transfusion set TS-203, Transfusion set TS-204
Product type: medical device
Design: sterile, ethylene oxide sterilization
Target patient: population of pediatric or adult patients
Designated user: medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used: PE, ABS, PVC, PS, PP, PA
Compatibility: With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With blood bags.
Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects: No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning: Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

The set is intended for therapeutic transferring of preserved blood, basic blood components, blood substitutes etc. from bottles or plastic bags into the patient.

Instruction for use

- Under aseptic conditions, remove the set from the packing.
- Close the hose regulation clamp.
- Take off the cover of the piercing spike and the cover of the cone.
- Inject the piercing spike upright into the plug of the transfusion bag. Hang the transfusion bag onto the stand.
- Pressing the drop chamber fill it by 1/3, open the regulation clamp, open the plug of the air inlet opening and let the air out from the hose.
- Close the regulation clamp.
- Connect the set on the patient, open the regulation clamp and start the transfusion.

Procedure for safe disposal of medical devices after use:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Filter of liquid with pore size TS-204
	The number of drops per milliliter		Medical Device		Date of Manufacture		Filter of liquid with pore size TS-201 and TS-203

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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