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**Product type:** medical device

**Design:** sterile, ethylene oxide sterilization

**Target patient:** population of pediatric or adult patients

**Designated user:** medical personnel with appropriate education, knowledge and practical experience in medical devices.

**Materials used:** PE, ABS, PVC, PS, PP, PA

**Compatibility:** With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With blood bags.

**Contraindication:** Do not use on any patient with a known allergic reaction to any of the product components.

**Adverse effects:** No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.

**Warning:** Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.


**Storage:** The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

**Intended use**

The set is intended for therapeutic transferring of preserved blood, basic blood components, blood substitutes etc. from bottles or plastic bags into the patient.

**Instruction for use**

1. Take the set out from the package.
2. Close the hose regulation clamp.
3. Take off the cover of the piercing spike and the cover of the cone.
4. Inject the piercing spike upright into the plug of the transfusion bag. Hang the transfusion bag onto the stand.
5. Pressing the drop chamber fill it by 1/3, open the regulation clamp, open the plug of the air inlet opening and let the air out from the hose.
6. Close the regulation clamp.
7. Connect the set on the patient, open the regulation clamp and start the transfusion.

**Procedure for safe disposal of medical devices after use:**

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

**Used symbols on packaging**

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use	Sterilised by Ethylene Oxide	
Pyrogen free		Keep away from sunlight	Latex free product	Does not contain DEHP phthalate	
Do not Resterilize		For single use only	Do not use if unit package is damaged	Filter of liquid with pore size TS-204	
The number of drops per milliliter		Medical Device	Date of Manufacture	Filter of liquid with pore size TS-201 and TS-203	

The instruction for use is also available in electronic form at [www.gama.cz](http://www.gama.cz). Free copies of other copies are available on request via e-mail: [sales@gama.cz](mailto:sales@gama.cz).

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