CS Trojcestný kohout M

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EN Three-way tap M

Catalogue No.:	Medical device variants:
V696434-M	Three-way tap – red NO PVC M
V696435-M	Three-way tap – blue NO PVC M
V696442-M	Three-way tap – white NO PVC M
Functionality:	the manufacturer declares the functionality of the device when used in accordance with the instructions for use
Duration of use:	„Short term“ - continuous use for the time maximum is 72 hours (3 days)
Design:	sterile, ethylene oxide sterilization
Target patient:	Patients undergoing intravenous or intra-arterial infusion therapy where concomitant administration of infusions or drugs is required.
Designated user:	Physicians, anesthetists, nurses, and other healthcare professionals in healthcare facilities who administer intravenous or intra-arterial infusion therapy.
Materials used:	PC, PE, PP, silicone fluid
Compatibility:	With medical devices equipped with standard Luer connectors under ČSN EN ISO 80369-7. With liquid drugs, infusion solutions or lipid emulsions. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects:	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning:	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient. If the tap is not exchanged after 72 hours of use, cracking may occur resulting in air embolism, leaks of the infusion solution and/or its contamination.
Storage and transport:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days

Intended use

Three-way taps are designed for concurrent infusion therapies when administering infusions or medications and can also be used for blood collection. They are intended for pressure and gravity applications up to a pressure of 2 bar for a maximum of 72 hours (3 days). They allow setting, controlling, controlling or stopping the desired direction of fluid flow. The three-way stopcock is compatible with common fluids used in infusion therapy, anesthesia, monitoring, as well as solutions that contain lipids. The color separation of the taps facilitates the identification of individual infusion lines or treatment applications. Three-way taps are compatible with Luer-Lock cones that meet the requirements of the ČSN EN ISO 80369-7 standard.

Instruction for use

- Under aseptic conditions remove the tap from the packaging and remove the guards.
- Attach the tap so as to ensure a safe connection.
- Turn the tap to set the desired direction of flow.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use	Sterilised by Ethylene Oxide and single sterile barrier system	
Pyrogen free		Keep away from sunlight	Latex free product	Ohne PVC	
Do not Resterilize		For single use only	Do not use if unit package is damaged	Pressure	
MD	Medical Device	Date and Country of Manufacture	UDI	Unique Device Identifier	

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