CS **Trojcestný kohout**

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Catalogue No.:	Medical device variants:
V696434	Three-way tap – red NO PVC
V696435	Three-way tap – blue NO PVC
V696442	Three-way tap – white NO PVC
Product type:	medical device
Design:	sterile, ethylene oxide sterilization
Target patient:	population of pediatric or adult patients
Designated user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used:	PC, PE
Compatibility:	With medical devices equipped with standard Luer connectors under ČSN EN ISO 80369-7. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects:	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning:	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient. If the tap is not exchanged after 72 hours of use, cracking may occur resulting in air embolism, leaks of the infusion solution and/or its contamination.
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

Three-way taps are designed for simultaneous infusion therapies to administer infusions or drugs and can also be used for blood collection. They are designed for gravity and pressure applications up to a pressure of 2 bar for up to 72 hours (3 days). They allow configuration, management, control or suspension of the fluid's desired flow direction. The three-way taps are compatible with all fluids used in infusion therapy, anaesthesia, monitoring as well as with solutions containing lipids. Colour coding of the taps makes it easy to identify individual infusion lines or therapeutic applications. The three-way taps are compatible with Luer-Lock cones, which meet the requirements of ČSN EN ISO 80369-7.

Instruction for use

- Under aseptic conditions remove the tap from the packaging and remove the guards.
- Attach the tap so as to ensure a safe connection.
- Turn the tap to set the desired direction of flow.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use	Sterilised by Ethylene Oxide	
Pyrogen free		Keep away from sunlight	Latex free product	Does not contain DEHP phthalate	
Do not Resterilize		For single use only	Do not use if unit package is damaged	Pressure	
MD	Medical Device	Date of Manufacture			

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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