Urodynamická souprava malá USM M

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EN Small urodynamic set USM M

Catalogue No.: V606481-02M
 V606481-03M
Medical device variants:
 Small urodynamic set USM -1845 mm M
 Small urodynamic set USM-S M
Functionality: the manufacturer declares the functionality of the device when used in accordance with the instructions for use
Duration of use: Short term" - continuous use for the time maximum is 24 hours
Design: sterile, ethylene oxide sterilization
Target patient: Patients with bladder dysfunction.
Designated user: Urologist or gynecologist.
Materials used: PVC, PE, ABS, PS, membrane VERSAPORE RC, PP, silicone
Compatibility: With urodynamic appliances UROMIC (including UROMIC SMABA), except for UROMIC 2 and UROMIC 4. With a cystometric urological catheter, which has a connection part consisting of a threaded end according to the ČSN EN ISO 80369-7 standard. With an infusion bag made of plastic, a plastic or glass bottle that has a rubber cap. With plastic infusion bags, plastic or glass bottles. Small urodynamic set USM is not intended for drug delivery.
Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects: No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning: Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage and transport: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action.
 The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!

Intended use

Single use medical device intended as accessory to UROMIC apparatus during urodynamic measurements. Small urodynamic set USM is designated for the UROMIC instrument series types, except for UROMIC 2 and UROMIC 4.

Use the sets only in connection with the prescribed type of pump.

Instruction for use

- Under aseptic conditions, remove the set from the packing. Close the pinchcock on the tube.
- Remove the cover from the puncturing style and the cover from the cone.
- Insert the puncturing style vertically into the stopper of the bottle or sac. Hand the bottle or sac on the stand.
- Fill up the drip chamber halfway.
- Open the stopper of the aeration opening, open also the pinchcock - and remove any air from the hose.

Note: If the set is used to transfer solution from the plastic pack, leave the aerating inlet closed.

- Close the pinchcock on the tube.
- Place the flexible part of the set into the pump, open the pinchcock.
- Connect the attachment cone to the opening of the catheter.
- Follow the instructions for using the pump for which the kit is designated.

NOTE: When the pump is running, the pinchcock must be fully open. After the dosing is finished, it is necessary to close the regulating pinchcock before the set is removed from the apparatus.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide and single sterile barrier system
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Pressure
	Medical Device		Date and Country of Manufacture		Unique Device Identifier		

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