## Urodynamická souprava malá USM

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## Small urodynamic set USM

<b>Catalogue No.:</b>	<b>Medical device variants:</b>
V606481-02ND	Small urodynamic set USM -1845 mm
V606481-03ND	Small urodynamic set USM-S
<b>Product type:</b>	medical device
<b>Design:</b>	sterile, ethylene oxide sterilization
<b>Target patient:</b>	population of pediatric or adult patients
<b>Designated user:</b>	medical personnel with appropriate education, knowledge and practical experience in medical devices.
<b>Materials used:</b>	PVC, PE, PP, PS, ABS, silicone
<b>Compatibility:</b>	With urodynamic appliances UROMIC (including UROMIC SMABA), except for UROMIC 2 and UROMIC 4. With cystometry urological catheter. With plastic infusion bags, plastic or glass bottles. Small urodynamic set USM is not intended for drug delivery.
<b>Contraindication:</b>	Do not use on any patient with a known allergic reaction to any of the product components.
<b>Adverse effects:</b>	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
<b>Warning:</b>	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
<b>Storage:</b>	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

### Intended use

Single use medical device intended as accessory to UROMIC apparatus during urodynamic measurements. Urodynamic kit small USM is designated for the UROMIC instrument series types, except for UROMIC 2 and UROMIC 4.  
**Use the sets only in connection with the prescribed type of pump.**

### Instruction for use

- Remove the set from the package. Close the pinchcock on the tube.
  - Remove the cover from the puncturing style and the cover from the cone.
  - Insert the puncturing style vertically into the stopper of the bottle or sac. Hand the bottle or sac on the stand.
  - Fill up the drip chamber halfway.
  - Open the stopper of the aeration opening, open also the pinchcock - and remove any air from the hose.
- Note: If the set is used to transfer solution from the plastic pack, leave the aerating inlet closed.
- Close the pinchcock on the tube.
  - Place the flexible part of the set into the pump, open the pinchcock.
  - Connect the attachment cone to the opening of the catheter.
  - Follow the instructions for using the pump for which the kit is designated.

**NOTE:** When the pump is running, the pinchcock must be fully open. After the dosing is finished, it is necessary to close the regulating pinchcock before the set is removed from the apparatus.

### PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

### Used symbols on packaging

REF	Catalogue Number		Keep away from Rain	LOT	Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Pressure
MD	Medical Device		Date of Manufacture				

The instruction for use is also available in electronic form at [www.gama.cz](http://www.gama.cz). Free copies of other copies are available on request via e-mail: [sales@gama.cz](mailto:sales@gama.cz).

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